



## Sartorius Microsart® e.motion Membrane Filter Dispensers SARTORIUS

Efficient and ergonomic workflows for reliable and fast membrane filter dispensing. Thanks to the small footprint the Microsart® e.motion dispenser saves precious space under the laminar air flow which can be used for samples, media and funnels instead. The Microsart® e.motion dispenser offers the easiest and most flexible way of membrane filter dispensing. The membrane filters are released from their sterile packaging fully automatically at the touch of a button or hands-free or by foot switch: the choice is yours!

- Capacity of up to 250 membrane filters
- Efficient workflow, fewer interruptions for reloading
- 3 operating modes
- Highest flexibility of operation
- Cableless operation and portability
- Easy to transport for convenient handling
- Dims, L x W x H: 233 x 172 x 230mm

Code	Alt Ref	Description	Price
<b>FIL1006</b>	16713----B0	With built in battery and AC adapter	<b>£1644.00</b>
<b>FIL1008</b>	16713----PS	With AC adapter	<b>£1358.00</b>



### Membrane Filters

Microsart® e.motion membrane filters give you the highest flexibility: they can be used with most commercially available membrane filter dispensers. Big packs containing 250 membrane filters enable you to work a complete shift without the need to stop and reload the dispenser.

- Sterile
- Pore size: 0.45µm high-flow
- Diameter: 47mm

**Supplied in packs of 1000 (4 x 250)**

Code	Alt Ref	Membrane Filter Colour	Grid Colour	Price
<b>FIL1012</b>	114H6Z-47---SJR	White	Black	<b>£499.00</b>
<b>FIL1014</b>	130H6Z-47---SJR	Grey*	White	<b>£499.00</b>

\* Grey membrane filters turn black when wet

## A world of products in one place...

- Over 500,000 products
- Easy ordering capability
- Detailed product information
- Promotions and offers

