



Whatman™ Mixed Cellulose Ester (WME) Membranes

cytiva
Whatman

Whatman mixed cellulose ester membranes are 47mm in diameter, and offer a variety of gridded formats for microbial colony quantification. Composed of cellulose acetate and cellulose nitrate, they offer a high degree of internal surface area for greater product adsorption.

Supplied in packs of 100

Code	Alt Ref	Description	Sterility	Pore Size, µm	Price
FIL3028	7153-104	Black, gridded, with pad	Sterile	0.45	£221.00
FIL3058	7187-114	White, gridded, without pad	Sterile	0.2	£97.00
FIL3070	7141-004	White, black grid	Non sterile	0.45	£158.00
FIL3072	7140-104	White, plain, with pad	Non sterile	0.45	£167.00
FIL3074	7141-104	White, gridded, with pad	Sterile	0.45	£95.00
FIL3076	7141-114	White, gridded, without pad	Sterile	0.45	£95.00
FIL3078	7141-204	White, black grid, autoclavable	Non sterile	0.45	£106.00



Whatman™ MicroPlus Cellulose Nitrate (CN) Membranes

cytiva
Whatman

Whatman offer a range of sterile MicroPlus cellulose nitrate (CN) membranes that are available folded in dispenser boxes compatible with most membrane dispensing systems (STL) or as single sterile packs (ST). STL membranes are widely used in the food, beverage, pharmaceutical and cosmetic industries and the membrane filtration method is excellent for the determination and enumeration of microbial contaminants.

- Developed for microbiological QC
- High mechanical strength
- High flow rate
- Pore size: 0.45µm
- Grid size: 3.1mm
- Type STL provided in dispenser boxes with 100 numbered membrane filters for easy removal and secure handling using Cytiva's Membrane-Butler (packs of 4 x 100)
- Type ST - single sterile packed (packs of 100)

Code	Alt Ref	Colour	Grid Colour	Ø, mm	Type	Pack	Price
FIL5420	10407112	White	Black	47	STL	4 x 100	£477.00
FIL5422	10407114	White	Black	50	STL	4 x 100	£498.00
FIL5424	10407713	White	Black	47	ST	1 x 100	£147.00
FIL5426	10407714	White	Black	50	ST	1 x 100	£151.00
FIL5428	10407132	Black	White	47	STL	4 x 100	£468.00
FIL5430	10407134	Black	White	50	STL	4 x 100	£498.00
FIL5432	10407734	Black	White	50	ST	1 x 100	£151.00
FIL5434	10407170	Green	Black	47	STL	4 x 100	£471.00
FIL5436	10407172	Green	Black	50	STL	4 x 100	£490.00



Whatman™ Membrane-Butler

cytiva
Whatman

The Membrane-Butler offers optimal handling of Whatman ME and MicroPlus CN STL membranes. With each turn a membrane filter is ejected from its sterile packaging and can be removed easily with tweezers.

- Durable design for reliable performance
- Easy grasp handle to dispense filter
- Minimised cross-contamination by eliminating handling of individually sealed filters
- Suitable for use on sterile benches
- Compact size for portability and preserve area of workspace

Code	Alt Ref	Price
FIL5475	10477100	£1430.00