



A&D GR Series Analytical Balances




Designed with a space saving small footprint. The unique door opening system allows both left or right handed operation via a single front mounted control making operation significantly quicker and easier for the operator.

- Auto self calibration
- 13 weighing units
- Data memory function
- One touch calibration
- Standard underhook
- Full GLP conformity (requires printer and interface)
- Easy access door system
- Splashproof keyboard and display
- Standard RS232C
- Dims, w x d x h: 249 x 330 x 327mm
- Pan dims: Ø 85mm

Code	Alt Ref	Capacity, g	Readability, mg	Price
BAL2060	GR-202-EC	42/210	0.01/0.1	£2421.00
BAL2064	GR-120-EC	120	0.1	£1648.00
BAL2062	GR-200-EC	210	0.1	£1700.00
BAL2066	GR-300	310	0.1	£2163.00

Accessories

Code	Alt Ref	Description	Price
BAL2800	AD-8127	Multifunctional compact printer	£649.00



IT'S IN OUR DNA...



0115 982 1111



scientificlabs.co.uk



**SCIENTIFIC
LABORATORY
SUPPLIES**