

Sartorius Microsart® Multi-Branch Manifolds Continued

SARTORIUS

Code	Alt Ref	Description	Price
FIL5232	168M1-MS	1 branch manifold for Microsart® disposables	£1627.00
FIL5234	168M2-MS	2 branch manifold for Microsart® disposables	£2490.00
FIL5236	168M3-MS	3 branch manifold for Microsart® disposables	£3023.00
FIL5238	168M6-MS	6 branch manifold for Microsart® disposables	£4623.00
FIL5240	168M1-BS100	1 branch manifold for Biosart® 100 monitors	£1652.00
FIL5242	168M2-BS100	2 branch manifold for Biosart® 100 monitors	£2515.00
FIL5244	168M3-BS100	3 branch manifold for Biosart® 100 monitors	£3074.00
FIL5246	168M6-BS100	6 branch manifold for Biosart® 100 monitors	£4712.00
FIL5248	168M1-BS250	1 branch manifold for Biosart® 250 monitors	£1546.00
FIL5254	168M3-BS250	2 branch manifold for Biosart® 250 monitors	£3023.00
FIL5256	168M6-BS250	3 branch manifold for Biosart® 250 monitors	£4674.00
FIL5258	168M1-SS100	6 branch manifold for Biosart® 250 monitors	£1716.00
FIL5260	168M3-SS100	3 branch manifold with 100mL stainless steel funnels	£3493.00
FIL5262	168M6-SS100	6 branch manifold with 100mL stainless steel funnels	£5664.00
FIL5264	168M1-SS500	1 branch manifold with 500mL stainless steel funnels	£1779.00
FIL5266	168M3-SS500	3 branch manifold with 500mL stainless steel funnels	£3671.00
FIL5268	168M6-SS500	6 branch manifold with 500mL stainless steel funnels	£5893.00

Whatman™ Vacuum Type Glass Membrane Holders



- Vacuum type glass holders produced from borosilicate glass
- Available with a choice of support screen
- Suitable for aqueous and organic solvent filtration
- Funnel seal ensures that the sample does not bypass the membrane and that particulates are retained on the surface of the membrane
- Sintered glass support is recommended for filtration and biological analysis
- 304 grade stainless steel support screen is suitable for use with proteinaceous solutions
- Comprises of reservoir, support screen, filter tube, stopper and clamp



Code	Alt Ref	Membrane Ø, mm	Filter Support	Capacity, mL	Price
FIL3200	1960-002	25	Borosilicate glass	25	£301.00
FIL3202	1960-032	25	Stainless steel	50	£403.00
FIL3204	1960-052	25	Stainless steel, mesh screen	25	£373.00
FIL3206	1960-004	47	Borosilicate glass	300	£396.00
FIL3208	1960-054	47	Stainless steel, mesh screen	300	£422.00

Visit our dedicated safety subsite

Bringing you easy to navigate categorisation, extensive choice and live pricing. Visit our safety subsite to view our ever increasing range of safety products.

scientificlabs.co.uk/safety (UK & Northern Ireland)

scientificlabs.ie/safety (Republic of Ireland)

scientificlabs.com/safety (International)

