



Whatman™ 3 Piece Filter Funnels

cytiva
Whatman

- Ideal for quick and easy filtration using high-efficiency glass microfibre filters
- Consists of a glass funnel and base, acrylic support platform, and 3 easy-to-apply metal clips.
- Range of diameters, funnel sizes, reservoir volumes and alternate support platforms provides a flexible system to meet the varied and unique needs of busy laboratories
- Chemically resistant construction of glass cylinder and base with a polymeric membrane support
- Provided with acrylic support and optional polypropylene or PTFE supports available separately for some sizes
- Clamps at edges to prevent peripheral loss due to fluid bypass around the edge of the filter
- Simplicity of design facilitates easy and thorough cleaning between use

| Code | Alt Ref | Reservoir Capacity, mL | Filter Paper Size, mm | Price |
|----------------|----------|------------------------|-----------------------|---------|
| FUN3100 | 1950-002 | 16 | 25 | £165.00 |
| FUN3102 | 1950-004 | 36 | 47 | £189.00 |
| FUN3104 | 1950-007 | 115 | 70 | £228.00 |
| FUN3106 | 1950-017 | 210 | 70 | £251.00 |
| FUN3108 | 1950-027 | 400 | 70 | £393.00 |
| FUN3110 | 1950-009 | 200 | 90 | £262.00 |
| FUN3112 | 1950-012 | 530 | 125 | £533.00 |

Alternative Support Plates

PTFE.

| Code | Alt Ref | Filter Paper Size, mm | Price |
|----------------|---------|-----------------------|--------|
| FUN3126 | 1950114 | 47 | £71.00 |
| FUN3130 | 1950117 | 70 | £90.00 |



SLS Select Syringe Filters



The highly asymmetrical pore structure of PES membranes offers an excellent loading capacity and high flow rate for our syringe filters which are naturally hydrophilic.

- Chemical stability: pH 4 – 8
- Recommended for: biology, preparation of protein samples, ultra-fast filtration of tissue culture media
- Diameter: 25mm
- Sterile

Supplied in packs of 50 (individually wrapped)

| Code | Pore Size, µm | Price |
|----------------|---------------|--------|
| SLS3011 | 0.22 | £67.00 |
| SLS3012 | 0.45 | £67.00 |