



SLS Lab Pro P-Series Precision Balances



- One touch automatic calibration
- Compact, B5 size footprint
- Fast, 1 second stabilisation
- Large bright vacuum fluorescent display
- Shock absorber function
- Multiple weighing units
- GLP/GMP/GCP/ISO compliant
- Statistical calculation function
- ACAI – Automatic Counting Accuracy Improvement
- Dims, w x d x h: 193 x 262.5 x 85.5mm

Code	Capacity, g	Readability, g	Pan Ø, mm	Price
SLS5540	122	0.001	130	£899.00
SLS5542	220	0.001	130	£927.00
SLS5062	320	0.001	130	£960.00
SLS5544	520	0.001	130	£1098.00
SLS5546	1220	0.01	150	£899.00
SLS5548	2200	0.01	150	£927.00
SLS5064	3200	0.01	150	£960.00
SLS5550	5200	0.01	150	£1098.00



SLS Lab Pro SK-i Series Compact Balances



- Single range, high resolution 1/6000 to 1/30000
- LCD display with backlight
- RS-232C as standard
- Auto power off function
- HI/OK/LO comparator indication
- Counting and percentage mode
- Multi unit: g, pieces, %, ounce, lb, oz, ct, mom, dwt, GN
- Optional Ni-MH rechargeable battery (metal hydrate)
- WinCT as standard
- ACAI counting function
- Dims, w x d x h: 190 x 218 x 551mm

Code	Capacity, g	Readability, g	Pan Dims, mm	Price
SLS5520	200	0.01	Ø 110	£282.00
SLS5522	300	0.01	Ø 110	£336.00
SLS5524	400	0.01	Ø 110	£4421.00
SLS5526	600	0.01	Ø 110	£512.00
SLS5528	1200	0.1	133 x 170	£265.00
SLS5530	2000	0.1	133 x 170	£282.00
SLS5532	3000	0.1	133 x 170	£347.00
SLS5534	4000	0.1	133 x 170	£431.00
SLS5535	6000	0.1	133 x 170	£473.00
SLS5536	6000	1.0	133 x 170	£282.00
SLS5538	12000	1.0	133 x 170	£392.00

Options and Accessories

Code	Description	Price
BAL2625	Rechargeable battery pack (Ni-MH)	£188.00
BAL2251	Carrying case	£87.00
BAL2800	Multifunctional compact printer	£649.00
BAL2572	External display	£244.00