



A

## SLS Select Conical Flasks



- Conical design allows for liquids to be stirred and swirled without being spilt
- Suitable for many different applications including titration
- With white graduations and marking spot
- Made from high quality borosilicate glass
- High chemical and thermal resistance
- Used for many laboratory applications
- Standards equivalent to ISO 1773

| Code           | Capacity, mL | Dims, Ø x h, mm | Neck Ø, mm | Pack | Price   |
|----------------|--------------|-----------------|------------|------|---------|
| <b>SLS5082</b> | 50           | 54 x 80         | 25         | 12   | £32.80  |
| <b>SLS5084</b> | 100          | 64 x 105        | 29         | 12   | £29.30  |
| <b>SLS5086</b> | 250          | 85 x 145        | 31         | 6    | £19.43  |
| <b>SLS5088</b> | 500          | 103 x 184       | 35         | 6    | £29.00  |
| <b>SLS5090</b> | 1000         | 131 x 214       | 42         | 6    | £41.40  |
| <b>SLS5286</b> | 3000         | 187 x 310       | 50         | 6    | £205.00 |
| <b>SLS5288</b> | 5000         | 220 x 365       | 50         | 2    | £102.00 |



A

## PYREX® Baffled Culture Flasks, Membrane Closure

**PYREX®**

- Ideal for aerobic microbiology and cell culture
- Complete with blue polypropylene GL45 screw cap with a gas-permeable PTFE membrane that permits gaseous exchange whilst contents remain sterile
- Can be autoclaved with the membrane screw cap fastened
- Manufactured from PYREX® borosilicate glass for maximum chemical and thermal resistance
- Four equidistant baffles provide consistent agitation when used with compatible orbital shakers
- Suitable for repeated autoclave cycles at 121°C

| Code           | Alt Ref | Capacity, mL | Height x Ø, mm | Pack | Price   |
|----------------|---------|--------------|----------------|------|---------|
| <b>FLA4302</b> | 1142/08 | 250          | 130 x 82       | 5    | £153.00 |
| <b>FLA4303</b> | 1142/10 | 500          | 165 x 101      | 5    | £161.00 |
| <b>FLA4304</b> | 1142/12 | 1000         | 205 x 129      | 2    | £141.00 |
| <b>FLA4305</b> | 1142/14 | 2000         | 265 x 166      | Each | £89.00  |

## Weigh up our range of balances!

Including portable, analytical and precision balances from OHAUS.

See pages 32–59

