

PYREX® Culture Flasks, Delong Neck


PYREX®

- Manufactured from PYREX® borosilicate glass for maximum chemical and thermal resistance
- Ideal for use in preparation of cultures
- Supplied without closures (available separately)
- Neck Ø: 38mm

Plain

- Compatible with shaker equipment
- White enamel markings

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Pack	Price
FLA4286	1135/14D	250	148 x 85	10	£195.00
FLA4288	1135/20D	500	183 x 105	5	£96.00
FLA4290	1135/26D	1000	229 x 131	2	£47.30
FLA4292	1135/30D	2000	302 x 166	Each	£31.50



A

Baffled

- Four equidistant baffles provide consistent agitation when used with compatible orbital shakers, causing a turbulent flow in the growth media increasing the oxygen intake of the cells
- Suitable for repeated autoclave cycles at 121°C

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Pack	Price
FLA4298	1141/08	250	155 x 82	5	£125.00
FLA4299	1141/10	500	193 x 101	5	£146.00
FLA4300	1141/12	1000	222 x 128	2	£110.00
FLA4301	1141/14	2000	300 x 166	Each	£79.00




A

Closures

- Two types of closure: anodised aluminium/2-position polypropylene
- Both closures are suitable for autoclaving
- Inner Ø: 38mm

Supplied in packs of 10

Code	Alt Ref	Material	Price
FLA4282	4550/38M	Aluminium	£110.00
FLA4284	4540/38P	Polypropylene	£33.50



A