



A

PYREX® Narrow Mouth Erlenmeyer Flasks

PYREX®

- Complies with ISO 1773
- Manufactured from PYREX® borosilicate glass for highest resistance to thermal shock and chemical corrosion
- White enamel graduations and marking spot
- Graduated to full capacity
- Printed with trace code for downloadable batch certificate

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA40600	1130/03	5*	52 x 25	12	10	£102.00
FLA40620	1130/05M	10*	60 x 30	12	10	£51.30
FLA40640	1130/11M	25	70 x 42	18	10	£46.10
FLA40660	1130/02M	50	85 x 51	18	10	£46.60
FLA40680	1130/06D	100	105 x 64	18	10	£48.40
FLA40670	1130/10D	150*	118 x 74	28	10	£56.00
FLA40690	1130/12D	200*	131 x 79	34	10	£56.00
FLA40700	1130/14D	250	145 x 85	34	10	£55.70
FLA40710	1130/16D	300*	156 x 87	34	10	£60.00
FLA40720	1130/20D	500	180 x 105	34	10	£70.00
FLA40740	1130/26D	1000	220 x 131	37	10	£119.00
FLA4076	1130/30D	2000	280 x 166	50	Each	£22.50
FLA4078	1130/32D	3000	310 x 187	52	Each	£41.20
FLA4080	1130/36D	5000	365 x 220	52	Each	£61.00

* 5, 10, 150, 200 & 300mL flasks are supplementary to the ISO range



A

PYREX® Narrow Mouth Erlenmeyer Flasks, Heavy Duty

PYREX®

- Mechanically stronger than standard flasks
- Designed for longer life and increased user safety
- Reinforced rims to reduce chipping and breakage

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA40840	4980-10	10	50 x 31	10	12	£169.00
FLA40860	4980-25	25	65 x 41	13	12	£120.00
FLA40880	4980-50	50	78 x 51	14	12	£120.00
FLA40900	4980-125	125	114 x 67	23	12	£102.00
FLA40920	4980-250	250	132 x 82	27	12	£110.00
FLA40940	4980-500	500	176 x 101	31	6	£62.00
FLA40960	4980-1L	1000	216 x 129	37	6	£100.00
FLA4098	4980-2L	2000	268 x 160	43	Each	£35.00
FLA4099	4980-4L	4000	360 x 206	43	Each	£95.00
FLA4097	4980-6L	6000	410 x 235	43	Each	£180.00