


PYREX® Wide Mouth Erlenmeyer Flasks

PYREX®

- Complies with ISO 24450
- Manufactured from PYREX® borosilicate glass
- Excellent thermal performance
- High resistance to chemical attack
- White enamel graduations and marking spot
- Printed with trace code for downloadable batch certificate

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA42500	1140/02M	50	85 x 48	34	10	£51.80
FLA42520	1140/04D	100	110 x 64	34	10	£53.20
FLA42540	1140/08D	250	140 x 85	50	10	£59.00
FLA42560	1140/10D	500	175 x 105	50	10	£79.00
FLA42580	1140/14D	1000	220 x 131	50	10	£141.00
FLA4260	1140/18D	2000	275 x 153	72	Each	£27.50



A

Simax Narrow Mouth Erlenmeyer Flasks


SIMAX®

- White enamel graduations and marking spot
- Borosilicate 3.3 glass

Plain

- Conform to ISO 1773

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA40000	24/25	25	75 x 42	22	10	£74.00
FLA40020	24/50	50	90 x 51	22	10	£46.40
FLA40040	24/100	100	105 x 64	22	10	£22.20
FLA40060	24/250	250	145 x 85	34	10	£24.60
FLA40080	24/500	500	180 x 105	34	10	£38.40
FLA40100	24/1000	1000	220 x 131	42	10	£58.00
FLA4012	24/2000	2000	280 x 166	50	Each	£11.70




A

Baffled

- Fluted indents down the length of the wall

Code	Alt Ref	Capacity, mL	Pack	Price
FLA41900	24/100/BAF	100	10	£186.00
FLA41920	24/250/BAF4	250	10	£196.00
FLA41940	24/500/BAF	500	5	£96.00
FLA4196	24/1000/BAF	1000	Each	£44.60
FLA4198	24/2000/BAF4	2000	Each	£61.00



A