


Simax Wide Mouth Erlenmeyer Flasks

SIMAX®

- White enamel graduations and marking spot
- Borosilicate 3.3 glass
- Conforms to ISO 24450

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA42000	25/50	50	85 x 51	34	10	£78.00
FLA42020	25/100	100	105 x 64	34	10	£25.20
FLA42040	25/250	250	140 x 85	50	10	£27.70
FLA42060	25/500	500	175 x 105	50	10	£42.00
FLA42080	25/1000	1000	220 x 131	50	10	£121.00
FLA4210	25/2000	2000	280 x 153	72	Each	£21.60




DURAN® Narrow Mouth Erlenmeyer Flasks

DURAN®

- ISO 1773:1997
- DURAN® borosilicate glass 3.3 to ISO 3585:1998; glass Class 1 to USP, EP & JP
- Tightly controlled glass quality, wall thickness and geometry ensure high resistance to thermal shock giving long service
- Very resistant to chemical attack
- Suitable for hotplate and autoclave use
- With marking label, graduations and retrace code which allows downloading of individual batch certificate (except 25mL)

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA40200	212161407	25	75 x 42	22	10	£55.00
FLA40220	212161707	50	90 x 51	22	10	£55.00
FLA40240	212162403	100	105 x 64	22	10	£55.00
FLA40260	212163202	200*	131 x 79	34	10	£55.00
FLA40280	212163605	250	145 x 85	34	10	£55.00
FLA40300	212163905	300*	156 x 87	34	10	£65.00
FLA40320	212164404	500	180 x 105	34	10	£74.00
FLA40340	212165409	1000	220 x 131	42	10	£124.00
FLA4036	212166302	2000	280 x 166	50	Each	£23.30
FLA4038	212166808	3000	310 x 187	52	Each	£42.50
FLA4040	212167307	5000	365 x 220	52	Each	£66.00

* 200 & 300mL additional to ISO 1773:1997



KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty