




A

DURAN® Round Flasks



- Flat bottom boiling flask to ISO 1773:1997
- Made from DURAN® borosilicate glass 3.3 to ISO 3585:1998; glass Class 1 to USP, EP & JP
- Tightly controlled glass quality, wall thickness and geometry ensure high resistance to thermal shock giving long service
- Very resistant to chemical attack
- Suitable for direct and indirect heating
- With marking label
- Wide mouth versions are available upon request

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA50200	217111708	50	90 x 51	26	10	£92.00
FLA50240	217113606	250	138 x 85	34	10	£82.00
FLA50260	217114405	500	163 x 105	34	10	£111.00
FLA50280	217115401	1000	190 x 131	42	10	£163.00
FLA5030	217116406	2000*	250 x 166	42	Each	£29.50
FLA5032	217116809	3000*	250 x 185	50	Each	£58.00
FLA5034	217117102	4000	275 x 207	50	Each	£83.00
FLA5036	217117308	5000*	290 x 223	50	Each	£94.00
FLA5038	217117608	6000	315 x 237	65	Each	£110.00
FLA5040	217118604	10000	360 x 280	65	Each	£199.00

* 2, 3 & 5L are additional to ISO 1773:1997



A

SLS Select Volumetric Flasks



- Manufactured to ISO 1042
- Borosilicate 3.3 glass
- Blue enamel graduations

Class A

- Supplied with PE plastic stopper and batch certificate

Code	Capacity, mL	Tolerance, ± mL	Stopper Size	Pack	Price
FLA77000	5	0.025	10/19	2	£8.10
FLA77020	10	0.025	10/19	2	£9.20
FLA77040	25	0.040	10/19	2	£9.70
FLA77060	50	0.060	12/21	2	£12.10
FLA77080	100	0.100	14/23	2	£22.40
FLA77100	200	0.150	14/23	2	£18.90
FLA77120	250	0.150	14/23	2	£19.80
FLA77140	500	0.250	19/26	2	£23.70
FLA77160	1000	0.400	24/29	2	£39.70
FLA7718	2000	0.600	29/32	Each	£31.00

Class B

- Supplied with PE plastic stopper

Code	Capacity, mL	Tolerance, ± mL	Stopper Size	Pack	Price
FLA70000	5	0.08	10/19	2	£35.00
FLA70020	10	0.08	10/19	2	£20.20
FLA70060	25	0.08	10/19	2	£9.30
FLA70080	50	0.12	12/21	2	£15.70
FLA70100	100	0.20	14/23	2	£12.60
FLA70120	200	0.3	14/23	2	£13.00
FLA70140	250	0.3	14/23	2	£17.40
FLA70160	500	0.5	19/26	2	£22.10
FLA70180	1000	0.8	24/29	2	£33.80
FLA7020	2000	2.4	29/32	Each	£27.20
FLA7022	5000	2.0	34/35	Each	£166.00