




A

PYREX® Amber Volumetric Flasks, Class A PYREX®

- Complies with ISO 1042 and DIN 12664
- Calibrated to the tolerance required by the United States Pharmacopoeia (USP) which are tighter than those of ISO/DIN standards
- Manufactured from borosilicate glass with an amber coating on exterior surfaces to facilitate use with light sensitive solutions
- Each flask is laser etched with an individual batch code and date of manufacture
- Conformity 'H' certified
- Batch certificate available to download

Code	Alt Ref	Capacity, mL	Tolerance, ± mL	Stopper Size	Pack	Price
FLA1140	3152/04A	5	0.020	10/19	5	£135.00
FLA1142	3152/06A	10	0.020	10/19	5	£135.00
FLA1144	3152/08A	25	0.03	10/19	5	£89.00
FLA1146	3152/10A	50	0.05	12/21	5	£101.00
FLA1148	3152/12A	100	0.08	14/23	5	£117.00
FLA1150	3152/14A	200	0.10	14/23	5	£158.00
FLA1152	3152/16A	250	0.12	14/23	2	£71.00
FLA1154	3152/18A	500	0.20	19/26	2	£92.00
FLA1156	3152/20A	1000	0.30	24/29	2	£127.00
FLA1158	3152/22A	2000	0.50	29/32	Each	£114.00




A

MBL® Volumetric Flasks, Class B MBL® VOLUMETRICS

- Complies with ISO 1042 and DIN 12664
- Manufactured from borosilicate glass for corrosion-free performance
- Inscriptions in blue enamel
- Polyethylene stopper

Code	Alt Ref	Capacity, mL	Tolerance, ± mL	Stopper Size	Pack	Price
FLA7230	FRD022P	5	0.05	10/19	5	£32.10
FLA7231	FRD024P	10	0.05	10/19	5	£32.10
FLA7232	FRD025P	20	0.08	10/19	5	£37.60
FLA7233	FRD026P	25	0.08	10/19	5	£37.70
FLA7234	FRD028P	50	0.12	12/21	5	£38.70
FLA7235	FRD030P	100	0.20	14/23	5	£47.70
FLA7236	FRD032P	200	0.30	14/23	5	£66.00
FLA7237	FRD034P	250	0.30	14/23	2	£26.60
FLA7238	FRD036P	500	0.50	19/26	2	£35.70
FLA7239	FRD038P	1000	0.80	24/29	2	£47.30
FLA7240	FRD040P	2000	1.20	29/32	Each	£41.30



A

Azlon® Volumetric Flasks, Class B AZLON®

- Translucent polypropylene with good chemical resistance
- Clear easy to read calibration line
- Calibrated to ISO 1042, class B tolerances
- Supplied with polypropylene stopper

Code	Alt Ref	Capacity, mL	Tolerance, ± mL	Stopper Size	Price
FLA7200	FDF012N	25	0.08	10/19	£14.00
FLA7202	FDF018N	50	0.12	10/19	£13.00
FLA7204	FDF024N	100	0.2	14/23	£14.00
FLA7206	FDF030N	250	0.3	19/26	£18.00
FLA7208	FDF036N	500	0.5	19/26	£20.50
FLA7207	FDF038N	1000	0.8	24/29	£22.00