

## Scanvac Laboratory Freeze Dryers Continued



### STEP 2 – Top Lid

Basics, Pro and Touch only.

Code	Alt Ref	Description	Price
<b>FRE4511</b>	7.001.100.061	ACpl manifold and chamber connection plate	<b>£374.00</b>



### STEP 3 – Manifold

Basics, Pro and Touch only.

Code	Description	Price
<b>FRE4512</b>	M4 basic manifold with 4 x 3/4" rubber valves	<b>£1480.00</b>
<b>FRE4513</b>	M4 extension manifold with 4 x 3/4" rubber valves	<b>£1109.00</b>
<b>FRE4560</b>	Drum manifold with 12 3/4" rubber valves (no ACpl lid manifold required)	<b>£2775.00</b>
<b>FRE4532</b>	AmpMan 16 ampoule manifold, 16 x 6mm pipes	<b>£458.00</b>

### STEP 4 – Chamber

Non-heated chamber (for use with the Basic, Pro & Touch), available in various sizes and formats. These include options for stoppering, attachment of flasks, and racks with removable shelves. Electrical heated chamber (for use with the Touch & Superior Touch XL & XS) that ensure uniform and faster drying, especially in the secondary drying stage. Giving increased recovery rates with controlled programmable energy input.

Note: extra racks, shelves and trays can also be supplied on request.



Code	Description	Heated	Price
<b>FRE4527</b>	Acrylic chamber with two adj. trays, 200mm Ø	Non-heated	<b>£972.00</b>
<b>FRE4528</b>	Acrylic chamber with three adj. trays, 300mm Ø	Non-heated	<b>£2220.00</b>
<b>FRE4529</b>	Stop-acrylic chamber with three shelves and stoppering, 300mm Ø	Non-heated	<b>£6408.00</b>
<b>FRE4530</b>	4V acrylic chamber with three adj. trays, 300mm Ø and 4 x 3/4" rubber valves	Non-heated	<b>£2775.00</b>
<b>FRE4562</b>	8V acrylic chamber with 3 stainless steel shelves, three adj. trays 300mm Ø and 8 x 3/4" rubber valves	Non-heated	<b>£3517.00</b>
<b>FRE4531</b>	Acrylic chamber with three adj. trays, 500mm Ø	Non-heated	<b>£7394.00</b>
<b>FRE4534</b>	Electrically heated acrylic chamber with five shelves, 300mm Ø	Electrically heated	<b>£8006.00</b>
<b>FRE4535</b>	Electrically heated acrylic chamber with five shelves and 4 x 3/4" rubber valves, 300mm Ø	Electrically heated	<b>£8559.00</b>
<b>FRE4536</b>	Electrically heated acrylic chamber with five shelves and 8 x 3/4" rubber valves, 300mm Ø	Electrically heated	<b>£9295.00</b>
<b>FRE4537</b>	Electrically heated acrylic chamber with three shelves and stoppering, 300mm Ø	Electrically heated	<b>£12241.00</b>
<b>FRE4564</b>	Electrically heated acrylic chamber with two shelves 500mm Ø	Electrically heated	<b>£13177.00</b>

### STEP 5 – Accessories

Code	Description	Price
<b>FRE4514</b>	C29/32 cone for freeze drying flasks	<b>£100.00</b>
<b>FRE4515</b>	C34/36 cone for freeze drying flasks	<b>£100.00</b>
<b>FRE4516</b>	C24/29 cone for freeze drying flasks	<b>£124.00</b>
<b>FRE4517</b>	RB 300mL freeze drying flask	<b>£235.00</b>
<b>FRE4518</b>	RB 500mL freeze drying flask	<b>£278.00</b>
<b>FRE4519</b>	RB 1000mL freeze drying flask	<b>£325.00</b>
<b>FRE4520</b>	RB 2000mL freeze drying flask	<b>£372.00</b>
<b>FRE4521</b>	150mL freeze drying flask/chamber with rubber lid and connector	<b>£223.00</b>
<b>FRE4522</b>	300mL freeze drying flask/chamber with rubber lid and connector	<b>£237.00</b>
<b>FRE4523</b>	600mL freeze drying flask/chamber with rubber lid and connector	<b>£350.00</b>
<b>FRE4524</b>	1200mL freeze drying flask/chamber with rubber lid and connector	<b>£463.00</b>

