

1
YEAR

Buchi Lyovapor™ L-200 Freeze Dryers

BUCHI

High-quality freeze drying (–55°C, 6kg) with a high level of automation. Infinite-Control™ includes easy method creation, method running, data logging, chart recording in real-time and interruption if needed from anywhere at anytime.

- Remote monitoring on mobile devices
- Stay informed with push notifications
- Save your precious product by the sample protective state
- User-friendly method handling and operation
- Various drying chamber options for your application needs
- Optional upgrade to Interface Pro available
- Easy to install on trolley, bench or in fume hood
- Optimise your process due to the end point determination
- Stable process parameters

| Code | Alt Ref | Description | Price |
|----------------|--------------------|---|-----------|
| FRE1054 | 11L200IN0102000001 | With 12 place manifold rack | £11256.00 |
| FRE1056 | 11L200IN0102001001 | With 12 place manifold rack and Pfeiffer duo pump | £13858.00 |
| FRE1058 | 11L200IN1412000001 | With 4 shelf chamber | £11384.00 |
| FRE1062 | 11L200IN1412001001 | With 4 shelf chamber and Pfeiffer duo pump | £13986.00 |

Accessories

| Code | Alt Ref | Description | Price |
|----------------|-------------|--|----------|
| FRE1064 | 11DCIN01000 | Drying rack 12 position manifold | £2314.00 |
| FRE1066 | 11DCIN02000 | Drying rack 24 position manifold | £3885.00 |
| FRE1068 | 11DCIN14100 | 4 PMMA shelf drying chamber | £2498.00 |
| FRE1072 | 11DCIN14200 | 4 PMMA shelf drying chamber and 12 position manifold | £2960.00 |
| FRE1074 | 11DCIN16100 | 6 PMMA shelf drying chamber | £3284.00 |
| FRE1076 | 11DCIN16200 | 6 PMMA shelf drying chamber and 12 position manifold | £3691.00 |
| FRE1078 | 11DCPH24110 | Heated 4 shelf drying chamber | £5941.00 |
| FRE1102 | 11DCPH24210 | Heated 4 shelf drying chamber and 12 position manifold | £6270.00 |
| FRE1106 | 11DCPH26110 | Heated 6 shelf drying chamber | £7180.00 |
| FRE1108 | 11DCPH26210 | Heated 6 shelf drying chamber and 12 position manifold | £7491.00 |
| FRE1112 | 11063692 | Stainless steel coated trolley | £1388.00 |



PYREX® Separating Funnels, Conical

PYREX®

- Complies to ISO 4800 type 1
- Manufactured from chemically resistant PYREX® borosilicate glass
- With interchangeable PTFE stopcock key
- Supplied with glass and polyethylene stopper

| Code | Alt Ref | Capacity, mL | Stopper Size | Stopcock Bore, mm | Price |
|----------------|----------|--------------|--------------|-------------------|---------|
| FUN4123 | 2190/01M | 50 | 19/26 | 1.5 | £70.00 |
| FUN4120 | 2190/02M | 100 | 19/26 | 1.5 | £70.00 |
| FUN4122 | 2190/04M | 250 | 24/29 | 2.5 | £77.00 |
| FUN4124 | 2190/06M | 500 | 24/29 | 4.0 | £81.00 |
| FUN4126 | 2190/08M | 1000 | 29/32 | 4.0 | £94.00 |
| FUN4128 | 2190/10M | 2000 | 29/32 | 6.0 | £126.00 |