



Büchner Funnels

- Porcelain
- Perforated support plate

Code	Capacity, mL	Filter Paper Size, mm	Top Ø, mm	Stem Length, mm	Price
FUN3002	70	55	55	64	£24.00
FUN3004	120	70	70	64	£29.00
FUN3006	240	90	90	71	£40.00
FUN3008	400	110	110	83	£52.30



Kartell Funnel Adapters

Kartell
LABWARE

- Supports for Büchner funnels
- Neoprene rubber

Includes the following:

- Ext Ø: 21 x 11mm, int Ø: 17 x 7mm, height: 21mm
- Ext Ø: 27 x 16mm, int Ø: 22 x 11mm, height: 22mm
- Ext Ø: 37 x 22mm, int Ø: 31 x 16mm, height: 25mm
- Ext Ø: 46 x 29mm, int Ø: 39 x 22mm, height: 29mm
- Ext Ø: 58 x 38mm, int Ø: 50 x 30mm, height: 35mm
- Ext Ø: 69 x 45mm, int Ø: 60 x 36mm, height: 40mm
- Ext Ø: 86 x 57mm, int Ø: 75 x 46mm, height: 45mm

Code	Alt Ref	Pack	Price
FUN3030	420	Set of 7	£51.00



Carbolite ELF Economy Furnaces

CARBOLITE
GERO 30-3000°C

The Carbolite Economy furnace has maximum temperature of 1100°C and incorporates the innovative design of the ELF range - a one-piece moulded low thermal mass ceramic fibre chamber and integral heating elements which allows the maximum temperature to be achieved in about half the time of a traditional "muffle" furnace design.

- NOT suited for ashing of organics or hydrocarbons (please see the Carbolite AAF series of ashing furnaces for this application)
- The Economy range incorporates a type 301 PID temperature with fully adjustable initial ramp heat up rate, and process timer
- Supplied with P301 controller, other controller options are available - please contact us for pricing and further information

Model	ELF 11/6	ELF 11/14	ELF 11/23
Max Temperature, °C	1100		
Capacity, L	6	14	23
Max Continuous Temperature, °C	1000	1000	1100
Internal Dims, h x w x d, mm	165 x 180 x 210	210 x 220 x 310	235 x 255 x 400
External Dims, h x w x d, mm (excluding chimney)	580 x 410 x 420	630 x 450 x 520	715 x 505 x 690
Maximum Power, W	2000	2600	5000
Holding Power, W	900	1300	1500
Nominal Heat Up Time (mins to 1000°C)	35	40	29
Temperature Sensor	Type K thermocouple		
Net Weight, kg	24	31	52

Code	Alt Ref	Model	Price
FUR1200	ELF1106-230SN-02-301	ELF 11/6	£1962.00
FUR1202	ELF1114-230SN-02-301	ELF 11/14	£2151.00
FUR1204	ELF1123-230SN-02-301	ELF 11/23	£3137.00