

Sartorius Entris® II Essential Line Precision Balances Continued

SARTORIUS

Code	Alt Ref	Capacity, g	Readability, mg	Pan Dims, mm	Windshield	Calibration	Price
BAL1178	BCE822-1S	820	10	182 x 182	—	External	£898.00
BAL1180	BCE1202-1S	1200	10	182 x 182	—	External	£990.00
BAL1182	BCE2202-1S	2200	10	182 x 182	—	External	£1090.00
BAL1184	BCE3202-1S	3200	10	182 x 182	—	External	£1190.00
BAL1186	BCE4202-1S	4200	10	182 x 182	—	External	£1280.00
BAL1188	BCE6202-1S	6200	10	182 x 182	—	External	£1470.00
BAL1190	BCE2201-1S	2200	100	182 x 182	—	External	£898.00
BAL1192	BCE5201-1S	5200	100	182 x 182	—	External	£990.00
BAL1194	BCE8201-1S	8200	100	182 x 182	—	External	£1090.00
BAL1196	BCE6200-1S	6200	1000	182 x 182	—	External	£898.00
BAL1198	BCE8200-1S	8200	1000	182 x 182	—	External	£990.00
BAL1002	BCE223i-1S	220	1	Ø 120	Manual	Internal	£1132.00
BAL1004	BCE323i-1S	320	1	Ø 120	Manual	Internal	£1232.00
BAL1006	BCE423i-1S	420	1	Ø 120	Manual	Internal	£1322.00
BAL1008	BCE623i-1S	620	1	Ø 120	Manual	Internal	£1512.00
BAL1024	BCE653i-1S	650	1	Ø 120	Top access	Internal	£1712.00
BAL1026	BCE622i-1S	620	10	182 x 182	—	Internal	£1042.00
BAL1028	BCE822i-1S	820	10	182 x 182	—	Internal	£1132.00
BAL1032	BCE1202i-1S	1200	10	182 x 182	—	Internal	£1232.00
BAL1034	BCE2202i-1S	2200	10	182 x 182	—	Internal	£1322.00
BAL1036	BCE3202i-1S	3200	10	182 x 182	—	Internal	£1625.00
BAL1038	BCE4202i-1S	4200	10	182 x 182	—	Internal	£1512.00
BAL1042	BCE6202i-1S	6200	10	182 x 182	—	Internal	£1712.00
BAL1044	BCE2201i-1S	2200	100	182 x 182	—	Internal	£1132.00
BAL1046	BCE5201i-1S	5200	100	182 x 182	—	Internal	£1232.00
BAL1048	BCE8201i-1S	8200	100	182 x 182	—	Internal	£1322.00
BAL1096	BCE6200i-1S	6200	1000	182 x 182	—	Internal	£1132.00
BAL1098	BCE8200i-1S	8200	1000	182 x 182	—	Internal	£1232.00

Sartorius Entris® II Advanced Line Precision Balances

SARTORIUS

No matter what you are weighing, the Entris® II Advanced balance range is always the right choice, offering unrivalled value and backed by almost 150 years of German engineering expertise. Featuring monolithic weigh cell technology for highly accurate results year in and year out, the fastest stabilisation time in its class using state of the art weighing sensors, a graphic touch display and 13 built in application programs.

- Built in real time level support
- Integrated protection systems with 3 configurable levels to determine valid weighing data and ensure only valid data is transferred to external devices
- High chemical resistance
- Hybrid screen for excellent readability and use combines intuitive, wear resistant LED and touch technology
- Controlled access to balance settings with user management
- Gapless documentation thanks to CalAuditTrail
- Enables multiple connection to external devices via two USB Type C interfaces
- Allows individualised printouts for up to 6 unique profiles
- Password protection for secure operation
- Verified balances are conformity assessed with EU type examination certificate



Code	Alt Ref	Capacity, g	Readability, mg	Pan Dims, mm	Windshield	Calibration	EU Verified	Price
BAL4576	BCA223-1S	220	1	Ø 120	Manual	External	No	£2245.00
BAL4588	BCA323-1S	320	1	Ø 120	Manual	External	No	£2423.00
BAL4598	BCA423-1S	420	1	Ø 120	Manual	External	No	£2634.00
BAL4618	BCA623-1S	620	1	Ø 120	Manual	External	No	£3045.00