



Gas Regulators

- Each unit has a cylinder contents gauge and a variable output pressure-indicating gauge
- Suitable for mounting on to standard UK pressurised gas cylinder
- Complete with 1/4" hoesetail

| Code | Description | Price |
|----------------|--------------------------------------|---------|
| GAS1035 | Carbon dioxide regulator, 0 - 4 bar | £172.00 |
| GAS1030 | Carbon dioxide regulator, 0 - 10 bar | £208.00 |
| GAS1040 | Hydrogen regulator 0 - 4 bar | £227.00 |
| GAS1041 | Hydrogen regulator 0 - 10 bar | £261.00 |
| GAS1050 | Nitrogen/air regulator 0 - 4 bar | £225.00 |
| GAS1043 | Nitrogen/air regulator 0 - 10 bar | £211.00 |
| GAS1060 | Oxygen regulator 0 - 4 bar | £208.00 |



Universal Gas Cylinder Trolley

- For the safe movement of 1, 2, 3 or 4 cylinders at one time (depending on cylinder size)
- This trolley is very manoeuvrable and makes light work of long distances
- The 3rd wheel is very quick and easy to deploy or stow
- Zinc plated steel frame
- 250mm MCP puncture proof rubber tyre main wheels with a 100mm elastic rubber 3rd wheel

| Code | Price |
|----------------|---------|
| GAS2050 | £450.00 |



Gas Cylinder Trolleys

- Manufactured specifically for oxygen cylinders in accordance with BS2718
- Cylinders sit on a sheet steel base
- Cushioned steel retaining straps with wing nuts secure the cylinders
- Anti-static wheels
- Front buffers protect adjacent equipment

| Code | For use with | Price |
|----------------|----------------------------|---------|
| GAS1200 | 1 x cylinder Ø up to 180mm | £135.00 |
| GAS1202 | 1 x cylinder Ø up to 300mm | £139.00 |



Gas Cylinder Bench Clamp

- Heavy-duty alloy cylinder holder, which can be secured to benches up to 70mm thick with two plated steel tommy screws
- Two countersunk holes are provided to accept wood screws for permanent fixing, complete with fire resistant belt and buckle
- Belt breaking strain: 750kg, stitching: 500kg
- Buckle is British Aerospace specification

| Code | For use with | Price |
|----------------|----------------------------|--------|
| GAS1110 | 1 x cylinder Ø up to 300mm | £82.00 |