

## Epredia™ CryoStar™ NX50 and NX70 Cryostats Continued

epredia | Authorized Distributor

Cryochamber Specifications	
Capacity, L	42
Chamber Temperature (@20°C), °C	-25 ± 2
Knife Holder Temperature (@20 to 30°C), °C	-27 to -20 ± 2
Active Object Head Cooling (@20°C)	-43 ± 2
Integrated Peltier Fast-Freezing Device (@20°C)	-57 ± 3

Code	Alt Ref	Description	Price
HIS3106	957220	CryoStar™ NX50 cryostat standard	£24242.00
HIS3108	957210	CryoStar™ NX50 cryostat with height adjustment	£25257.00
HIS3110	957140	CryoStar™ NX50 cryostat with cold disinfection	£26271.00
HIS3112	957200	CryoStar™ NX50 cryostat with vacutome	£26472.00
HIS3114	957160	CryoStar™ NX50 cryostat with vacutome and cold disinfection	£29314.00
HIS3116	957040	CryoStar™ NX70 cryostat	£37083.00
HIS3118	957050	CryoStar™ NX70 cryostat with cold disinfection	£38905.00
HIS3120	957060	CryoStar™ NX70 cryostat with vacutome	£41638.00
HIS3122	957070	CryoStar™ NX70 cryostat with vacutome and cold disinfection	£43460.00

## Epredia™ High Capacity Section Dryer

epredia | Authorized Distributor

Provide maximum drying capacity for high workload laboratories. With space for 13 slide racks (of 30 capacity), it provides ample space, but is still compact enough to be located in close proximity to the microtome.

- Small footprint
- Anti-microbial paint finish (silver ion) to inhibit growth of bacteria
- Fan-assisted heating reduces section drying time
- Ultra-fast heating system for rapid warm up
- Digital temperature control guarantees heat accuracy within the chamber

Code	Alt Ref	Temp Range, °C	Price
HIS3046	A84600050	Room temp: 85 (±1) at 25 ambient	£1962.00



## Epredia™ Digital Round Tissue Section Bath

epredia | Authorized Distributor

Features innovative digital controls that maintain water temperatures, allowing tissue sections to be efficiently flattened prior to being lifted onto a microscope slide. It has an extra deep bowl, allowing slides to be fully dipped at any angle without hitting the bottom.

- Deep bowl for easier section pick up
- Wide bowl accommodates longer section ribbons
- Large area to rest hands when floating sections
- Miniature circuit breakers for ease-of-use; simply push to reset
- Anti-microbial paint finish (silver ion) to inhibit bacterial growth
- Ultra-fast heating system
- Significant slide drying area at the rear of the instrument
- A secure-fitting lid reduces airborne contamination, evaporation, and improves heating efficiency
- Automatic over-temperature cut-out for safety

Code	Alt Ref	Capacity, L	Price
HIS3044	A84600060	3.5	£1163.00

