

## Simport® Scientific Embedding Cassettes



Disposable plastic tissue cassettes suitable for holding and identifying tissue samples in processing, embedding and sectioning procedures.


- Fit securely in microtome chuck adapters
- Resistant to the chemical action of most histological solvents
- Slanted writing surface accepts markings easily, permitting sample identification throughout all stages of embedding and long afterwards when in archives
- Made of acetal
- Metal lid sold separately (HIS0020)
- Other colours are available, please enquire

Supplied in packs of 1500 (3 dispenser boxes of 500 cassettes)

Code	Alt Ref	Colour	Price
HIS0009	M480-2	White	£175.00
HIS0010	M480-3	Pink	£175.00
HIS0011	M480-4	Green	£175.00
HIS0012	M480-5	Yellow	£175.00

### Metal Lids

Code	Alt Ref	Pack	Price
HIS0020	M481	25	£82.00



## Simport® Scientific Histosette® I Cassettes



Disposable plastic cassettes specially designed to hold biopsy or tissue specimens during the embedding process as well as in a storage file. Moulded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents.

- One piece integral lid eliminates the need for separate steel lids
- Just snap apart and conveniently lock the lid into the base of the cassette
- Can be opened or closed as often as necessary and they always relock securely without danger of specimen loss
- Anterior writing area at a 30° angle
- Made of acetal plastic

Supplied in cases of 1500 (3 dispenser boxes of 500 cassettes)

Code	Alt Ref	Description	Colour	Price
HIS0021	M490-2	Tissue cassettes	White	£160.00
HIS0022	M490-3	Tissue cassettes	Pink	£160.00
HIS0023	M490-4	Tissue cassettes	Green	£160.00
HIS0024	M490-5	Tissue cassettes	Yellow	£160.00
HIS0025	M490-6	Tissue cassettes	Blue	£160.00
HIS0026	M490-7	Tissue cassettes	Peach	£160.00
HIS0027	M490-8	Tissue cassettes	Tan	£160.00
HIS0028	M490-9	Tissue cassettes	Grey	£160.00
HIS0029	M490-10	Tissue cassettes	Lilac	£160.00
HIS0030	M490-11	Tissue cassettes	Orange	£160.00
HIS0031	M490-12	Tissue cassettes	Aqua	£160.00
HIS0032	M491-2	Biopsy cassettes	White	£160.00
HIS0033	M491-3	Biopsy cassettes	Pink	£160.00
HIS0034	M491-4	Biopsy cassettes	Green	£160.00
HIS0035	M491-5	Biopsy cassettes	Yellow	£160.00
HIS0036	M491-6	Biopsy cassettes	Blue	£160.00
HIS0037	M491-7	Biopsy cassettes	Peach	£160.00
HIS0038	M491-8	Biopsy cassettes	Tan	£160.00
HIS0039	M491-9	Biopsy cassettes	Grey	£160.00
HIS0040	M491-10	Biopsy cassettes	Lilac	£160.00
HIS0041	M491-11	Biopsy cassettes	Orange	£160.00
HIS0042	M491-12	Biopsy cassettes	Aqua	£160.00

