



Simport® Scientific Macrosette™ Cassettes



Disposable plastic cassettes designed to hold larger tissue specimens during the embedding process, as well as in a storage cabinet. Same footprint as the Histosette® cassette but the Macrosette® is twice as high.

- One-piece integral lid eliminates the need for separate steel lids
- Made of acetal
- Specially made to hold larger tissue specimens
- Fully usable anterior writing area providing more labelling space
- Easy one hand manipulation

Supplied in case of 750 (3 dispenser boxes of 250 cassettes)

Code	Alt Ref	Dims, w x d x l, mm	Colour	Price
HIS0219	M512	28.5 x 13 x 40.1	White	£300.00



Simport® Scientific Stainless Steel Base Moulds



Stainless steel base moulds fit most histology cassettes. Superior thermal exchange with rounded corners. These moulds are for all applications in specimen embedding with Simport embedding rings and cassettes.

- Highly polished surface for easy paraffin block removal
- Rounded well corners for optimal paraffin ribboning
- Optimal thermal conductivity

Supplied in packs of 12

Code	Alt Ref	Size, mm	Price
HIS2008	M474-1	7 x 7 x 5	£102.00
HIS2010	M474-2	15 x 15 x 5	£102.00
HIS2012	M474-3	24 x 24 x 5	£102.00
HIS2014	M474-4	30 x 24 x 5	£102.00
HIS2016	M474-5	37 x 24 x 5	£102.00



Simport® Scientific Deep Base Moulds



Disposable base moulds offer ease and convenience of operation. They are inexpensive enough to be discarded after use, yet strong enough to be reused. Designed especially to be compatible with the HIS0219 (M512) Macrosette®.

- Smooth interior finish
- Rounded corners facilitating specimen removal
- Excellent thermal exchange
- Disposable
- Made of PVC

Code	Alt Ref	Dims, w x d x h, mm	Volume, mL	Pack	Price
HIS0225	M475-10	37 x 24 x 10	5	500	£102.00