

## OHAUS PR Precision Balances Continued



| Code    | Alt Ref  | Capacity, g | Readability, g | Pan Dims, mm | Windshield | Calibration | Price   |
|---------|----------|-------------|----------------|--------------|------------|-------------|---------|
| BAL1488 | PR223/E  | 220         | 0.001          | Ø 120        | Manual     | External    | £625.00 |
| BAL1490 | PR223    | 220         | 0.001          | Ø 120        | Manual     | Internal    | £765.00 |
| BAL1492 | PR423/E  | 420         | 0.001          | Ø 120        | Manual     | External    | £640.00 |
| BAL1494 | PR423    | 420         | 0.001          | Ø 120        | Manual     | Internal    | £780.00 |
| BAL1496 | PR523/E  | 520         | 0.001          | Ø 120        | Manual     | External    | £675.00 |
| BAL1498 | PR523    | 520         | 0.001          | Ø 120        | Manual     | Internal    | £820.00 |
| BAL1548 | PR1602/E | 1600        | 0.01           | Ø 180        | —          | External    | £640.00 |
| BAL1562 | PR1602   | 1600        | 0.01           | Ø 180        | —          | Internal    | £780.00 |
| BAL1564 | PR2202/E | 2200        | 0.01           | Ø 180        | —          | External    | £660.00 |
| BAL1566 | PR2202   | 2200        | 0.01           | Ø 180        | —          | Internal    | £810.00 |
| BAL1568 | PR4202/E | 4200        | 0.01           | Ø 180        | —          | External    | £690.00 |
| BAL1570 | PR4202   | 4200        | 0.01           | Ø 180        | —          | Internal    | £840.00 |
| BAL1572 | PR4201/E | 4200        | 0.1            | Ø 180        | —          | External    | £585.00 |
| BAL1574 | PR4201   | 4200        | 0.1            | Ø 180        | —          | Internal    | £730.00 |
| BAL1576 | PR6201/E | 6200        | 0.1            | Ø 180        | —          | External    | £610.00 |
| BAL1578 | PR6201   | 6200        | 0.1            | Ø 180        | —          | Internal    | £750.00 |



## OHAUS Scout® STX Precision Balances



Ready for your lab, wherever that may be. Packed with new design features including a graphic LCD touchscreen, the OHAUS Scout® Precision balances give exceptional weighing results and the choice of either mains or battery operation.

- Applications include weighing (7 units), counting, percent, display hold, totalisation, check weigh, specific gravity/density, dynamic
- 4.3" full colour TFT graphic LCD touchscreen
- Advanced weighing technology allows for one-second stabilisation time for increased productivity
- Integrated overload protection system protects the balance in demanding environments
- Mains or battery
- ABS plastic housing with stainless steel pan

| Code    | Alt Ref | Capacity, g | Readability, g | Pan Dims, mm | Windshield | Price   |
|---------|---------|-------------|----------------|--------------|------------|---------|
| BAL0041 | STX123  | 120         | 0.001          | Ø 93.2       | Hinged     | £585.00 |
| BAL0042 | STX223  | 220         | 0.001          | Ø 93.2       | Hinged     | £660.00 |
| BAL0043 | STX222  | 220         | 0.01           | Ø 120        | —          | £418.00 |
| BAL0044 | STX422  | 420         | 0.01           | Ø 120        | —          | £515.00 |
| BAL0045 | STX622  | 620         | 0.01           | Ø 120        | —          | £555.00 |
| BAL0046 | STX1202 | 1200        | 0.01           | 170 x 140    | —          | £620.00 |
| BAL0047 | STX2202 | 2200        | 0.01           | 170 x 140    | —          | £660.00 |
| BAL0048 | STX421  | 420         | 0.1            | Ø 120        | —          | £256.00 |
| BAL0049 | STX621  | 620         | 0.1            | 170 x 140    | —          | £293.00 |
| BAL0050 | STX2201 | 2200        | 0.1            | 170 x 140    | —          | £418.00 |
| BAL0051 | STX6201 | 6200        | 0.1            | 170 x 140    | —          | £560.00 |
| BAL0052 | STX8200 | 8200        | 1              | 170 x 140    | —          | £429.00 |

## Accessories

| Code    | Alt Ref  | Description             | Pack | Price   |
|---------|----------|-------------------------|------|---------|
| BAL1018 | 30268985 | Bluetooth interface kit | Each | £104.50 |
| BAL1019 | 30268987 | Stacking kit            | 6    | £50.10  |
| BAL1017 | 30268982 | RS232 interface kit     | Each | £55.70  |