

SLS Select Tissue Grinders



Supplied fully assembled, gamma sterilised and in easy-to-open Tyvek® package. Supplied with an extra cap for storage and transport, plus a patient identification label for traceability.




Abrasive Tip Grinders

Engineered for safety and ease of use with a patented closed system design that provides the ability to grind the specimen in a sealed container helping to prevent any risk of contact with the sample. The narrow pestle eliminates the problem of sample overflow common with other tissue grinders. The moulded-in abrasive surface does not leave any sediment to obstruct the sample during examination.

Supplied in packs of 10

Code	Capacity, mL	Description	Price
SLS1350	15	Conical	£131.00
SLS1352	35	Free-standing	£131.00
SLS1354	50	Conical	£131.00



Ultra Grinders

The pestle tips have been vitrified using a proprietary silica casting process, resulting in glass-like abrasive surface which provides a more aggressive grinding action to enable efficient processing of difficult and fibrous samples. Features the same patented closed system as standard grinders.

Supplied in packs of 10

Code	Capacity, mL	Description	Price
SLS1360	15	Conical	£150.00
SLS1362	50	Conical	£150.00



Griffiths Homogeniser Tissue Grinders

Produced from standard tubing unlike the majority of tissue grinders which are produced from precision bore tubing. These homogenisers are perfect for general work where tight tolerances are unnecessary.

Code	Description	Capacity, mL	Chamber Length, mm	Overall Length, mm	Pestle Ø, mm	Price
HOM3560	For grinding bacterial emulsions and tissue	5	80	150	12	£33.40
HOM3562	For coarse grinding	15	95	305	15	£47.20