

IKA ULTRA-TURRAX® Homogenisers Continued



T 25 Dispersing Tools

Code	Alt Ref	Description	Material	Max Temp, °C	Working Volume, mL	Stator/Rotor Ø, mm	Pack	Price
HOM2100	0001024200	S 25 N - 8 G	PTFE, AISI 316L	180	1 - 50	8/6.1	Each	£1717.00
HOM2492	0004446900	S 25 N - 8 G ST, saw tooth	PTFE, AISI 316L	180	1 - 50	8/6.1	Each	£1616.00
HOM2442	0000594000	S 25 N - 10 G	PTFE, AISI 316L	180	1 - 100	10/7.5	Each	£1287.00
HOM2494	0004447100	S 25 N - 10 G ST, saw tooth	PTFE, AISI 316L	180	1 - 100	10/7.5	Each	£1256.00
HOM2464	0003452200	S 25 D - 10 G - KS	Polycarbonate, polyetheretherketon	100	10 - 100	10/6.75	12	£362.00
HOM2466	0003452100	S 25 D - 14 G - KS	Polycarbonate, polyetheretherketon	100	10 - 500	14/9.5	12	£362.00
HOM2074	0020002973	S 25 KR - 18 G ST, saw tooth	AISI 316L, PTFE, FFKM	150	10 - 1500	18/13.4	Each	£1616.00
HOM2072	0020002971	S 25 KR - 18 G	AISI 316L, PTFE, FFKM	150	10 - 1500	18/12.7	Each	£1674.00
HOM2602	0002348000	S 25 KV - 18 G	FFPM, SIC, AISI 316L	220	10 - 1500	18/12.7	Each	£2446.00
HOM2444	0000593400	S 25 N - 18 G	PTFE, AISI 316L	180	10 - 1500	18/12.7	Each	£1178.00
HOM2496	0004447300	S 25 N - 18 G ST, saw tooth	PTFE, AISI 316L	180	10 - 1500	18/13.4	Each	£1263.00
HOM2075	0020002974	S 25 KR - 25 G ST, saw tooth	AISI 316L, PTFE, FFKM	150	50 - 2000	25/20	Each	£1929.00
HOM2073	0020002972	S 25 KR - 25 G	AISI 316L, PTFE, FFKM	150	50 - 2000	25/17	Each	£2004.00
HOM2606	0002466900	S 25 KV - 25 G	FFPM, SIC, AISI 316L	220	50 - 2000	25/17	Each	£2575.00
HOM2106	0001713300	S 25 N - 25 G	PTFE, AISI 316L	180	50 - 2000	25/17	Each	£1645.00
HOM2498	0004447500	S 25 N - 25 G ST, saw tooth	PTFE, AISI 316L	180	50 - 2000	25/20	Each	£1616.00
HOM2076	0020002975	S 25 KR - 25 F	AISI 316L, PTFE, FFKM	150	100 - 2000	25/18	Each	£1973.00
HOM2604	0002404000	S 25 KV - 25 F	FFPM, SIC, AISI 316L	220	100 - 2000	25/18	Each	£2514.00
HOM2108	0001713800	S 25 N - 25 F	PTFE, AISI 316L	180	100 - 2000	25	Each	£1616.00

T 50 Dispersing Tools

Code	Alt Ref	Description	Material	Max Temp, °C	Working Volume, mL	Stator/Rotor Ø, mm	Pack	Price
HOM2080	0008004100	S 50 KG - HH - G 45 F	FFPM, WC, AISI 316L	180	250 - 10000	45/40	Each	£5603.00
HOM2071	0010000631	S 50 KR - G 45 F	FFPM, PTFE, AISI 316L	150	250 - 10000	45/40	Each	£5030.00
HOM2448	0008003900	S 50 N - G 45 F	PTFE, AISI 316L	180	250 - 10000	45/40	Each	£3875.00
HOM2430	0001689300	R 50 high speed stirring shaft	Stainless steel (AISI 316L)	—	250 - 30000	—	Each	£2447.00
HOM2070	0010000630	S 50 KR - G 45 M	FFPM, PTFE, AISI 316L	150	500 - 15000	45/40.5	Each	£4845.00
HOM2079	0008003500	S 50 KG - HH - G 45 M	FFPM, WC, AISI 316L	180	500 - 15000	45/40.5	Each	£5421.00
HOM2446	0008003300	S 50 N - G 45 M	PTFE, AISI 316L	180	500 - 15000	45/40.5	Each	£3695.00
HOM2077	0008003200	S 50 KG - HH - G 45 G	FFPM, WC, AISI 316L	180	500 - 20000	45/36	Each	£5082.00
HOM2078	0010000648	S 50 KG - HH - G 45 G ST, saw tooth	FFPM, SIC, AISI 316L	180	500 - 20000	45/36	Each	£5082.00
HOM2068	0010000629	S 50 KR - G 45 G	FFPM, PTFE, AISI 316L	150	500 - 20000	45/36	Each	£4510.00
HOM2069	0010000632	S 50 KR - G 45 G ST, saw tooth	FFPM, PTFE, AISI 316L	150	500 - 20000	45/38	Each	£4510.00
HOM2202	0008003000	S 50 N - G 45 G	PTFE, AISI 316L	180	500 - 20000	45/36	Each	£3358.00
HOM2504	0008039500	S 50 N - G 45 G ST, saw tooth	PTFE, AISI 316L	180	500 - 20000	45/38	Each	£3358.00
HOM2610	0008005100	S 50 N - W 65 SK, cutting head	PTFE, AISI 316L	180	1000 - 10000	65/55	Each	£3765.00
HOM2612	0008006300	S 50 N - W 80 SMK, jet mixer head	PTFE, AISI 316L	180	1000 - 50000	80/34	Each	£3523.00