

OHAUS Explorer® Precision Balances



There are four user programmable touch free sensors which help to minimise the risk of cross contamination and avoid the need to repeatedly touch the balance.

- High resolution multilingual colour VGA touch display with icons for simple operation
- Programmable multifunctional touch free sensors for cleanliness, speed and safety
- Sturdy modular design includes metal housing and a display that can be separated from the weighing base and tower or wall mounted for easy viewing in any workspace
- Over 14 application modes including pipette adjustment
- Library for storage of personal application settings
- RS232 and USB as standard, plus 3rd connection slot (ethernet or 2nd RS232)
- Full GLP/GMP including date and time
- Supplied with data transfer software as standard
- Modular display design with adjustable angle (attach to base, wall or optional tower)
- Fast stabilisation times, optimised linearity and repeatability specifications, and enhanced vibration filtering
- Available with or without AutoCal™ automatic calibration system
- Level assist feature
- Integral weigh below hook
- All 1mg models come with 3 door draftshield



Code	Alt Ref	Model	Capacity, g	Readability, mg	Pan Dims, mm	Calibration	Price
BAL5111	30061979	EX223/E	220	1	Ø 130	External	£2300.00
BAL5110	83021341	EX223	220	1	Ø 130	Internal	£2600.00
BAL5113	30061980	EX423/E	420	1	Ø 130	External	£2905.00
BAL5112	83021342	EX423	420	1	Ø 130	Internal	£3205.00
BAL5116	83021348	EX623	620	1	Ø 130	Internal	£3660.00
BAL5118	83021350	EX1103	1100	1	Ø 130	Internal	£4610.00
BAL5123	30061981	EX2202/E	2200	10	180 x 200	External	£2200.00
BAL5122	83021356	EX2202	2200	10	180 x 200	Internal	£2505.00
BAL5125	30061982	EX4202/E	4200	10	180 x 200	External	£2850.00
BAL5124	83021357	EX4202	4200	10	180 x 200	Internal	£3155.00
BAL5135	30061983	EX6202/E	6200	10	180 x 200	External	£3370.00
BAL5128	83021360	EX6202	6200	10	180 x 200	Internal	£3665.00
BAL5130	83021362	EX10202	10200	10	180 x 200	Internal	£4310.00
BAL5133	30061984	EX6201/E	6200	100	180 x 200	External	£2050.00
BAL5134	83021361	EX6201	6200	100	180 x 200	Internal	£2345.00
BAL5136	83021365	EX10201	10200	100	180 x 200	Internal	£2885.00



OHAUS Adventurer® Precision Balances



Striking the ideal balance between inventive features and functional, uncomplicated weighing capabilities, the OHAUS Adventurer® incorporates all of the applications necessary for routine weighing and measurement activities. With a colour touchscreen, GLP/GMP compliance capabilities, and two USB ports.

- Wide viewing angle full colour touch screen
- Dual USB ports with GLP/GMP tracking capabilities
- Available with or without AutoCal™ automatic calibration system
- 1mg models have a two piece, top mounted draftshield with doors that open wide, providing ample room for sample placement