

IKA ULTRA-TURRAX® Homogenisers Continued



T 65 Dispersing Tools

| Code | Alt Ref | Description | Material | Max Temp, °C | Working Volume, mL | Stator/Rotor Ø, mm | Pack | Price |
|----------------|------------|-----------------------|-----------------------|--------------|--------------------|--------------------|------|-----------------|
| HOM2454 | 0008005900 | S 65 KG - HH - G 65 F | FFPM, WC, AISI 316L | 180 | 2000 - 30000 | 65/58 | Each | £5794.00 |
| HOM2087 | 0010000643 | S 65 KR - G 65 F | FFPM, PTFE, AISI 316L | 180 | 2000 - 30000 | 65/58 | Each | £5262.00 |
| HOM2452 | 0008005700 | S 65 KG - HH - G 65 M | FFPM, WC, AISI 316L | 180 | 2000 - 40000 | 65/58 | Each | £5577.00 |
| HOM2086 | 0010000642 | S 65 KR - G 65 M | FFPM, PTFE, AISI 316L | 180 | 2000 - 40000 | 65/58 | Each | £5048.00 |
| HOM2450 | 0008005500 | S 65 KG - HH - G 65 G | FFPM, WC, AISI 316L | 180 | 2000 - 50000 | 65/53 | Each | £5226.00 |
| HOM2085 | 0010000641 | S 65 KR - G 65 G | FFPM, PTFE, AISI 316L | 180 | 2000 - 50000 | 65/53 | Each | £4695.00 |



Cole-Parmer® SHM3 Digital Homogeniser

Cole-Parmer®

Can be stand mounted for larger samples, or used as a handheld unit for processing smaller batches. Digital speed control helps increase repeatability, and the soft start motor reduces splashing, keeping both you and your samples safer. Utilise the built in timer to process for seconds to minutes – also features a 30 second countdown for remote applications.

- Powerful and versatile for all your homogenisation needs
- Ideal for tissue homogenisation, cell disruption, emulsions, suspensions, and protein extraction
- Digital control for better repeatability and accuracy
- Process samples from 200µL to 2L
- Vary the 850W motor from 5000 to 25000RPM

| Code | Alt Ref | Model | Sample Size, mL | Speed Range, RPM | Dims, h x Ø, mm | Weight, kg | Price |
|----------------|----------|-------|-----------------|------------------|-----------------|------------|-----------------|
| HOM2508 | 04720-00 | SHM3 | 0.2 - 2000 | 5000 - 25000 | 70 x 216 | 1.6 | £2199.00 |

Rotor Probe

| Code | Alt Ref | Capacity, mL | Dims, l x Ø, mm | Price |
|----------------|----------|--------------|-----------------|----------------|
| HOM2510 | 04720-03 | 0.03 - 5 | 75 x 5 | £919.00 |

KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty