



Benchmark Digital Hotplates



Reliable performance and safe operation for everyday laboratory use. Featuring a durable, chemical resistant white ceramic top plate and advanced microprocessor controls, the constant monitoring system provides precise control of temperature with rapid heating and faster boiling times.

- Temperature range: ambient +5 to 400°C
- Temperature accuracy: $\pm 0.3^\circ\text{C}$
- Temperature viewed on the large digital display
- Safety LED illuminates when the top plate temperature exceeds 50°C
- Large, ceramic work surface
- Control actual sample temperature (with optional probe – sold separately)

Code	Alt Ref	Dims, w x d x h, mm	Platform Dims, mm	Weight, kg	Price
MIX1262	H3770-H-E	180 x 279 x 101	178 x 178	2.9	£291.00
MIX1270	H3710-H-E	295 x 424 x 114	254 x 254	7.5	£459.00



Mettler Toledo Densito Handheld Density Meters



Made for use in the lab and on the go, the Densito is perfect for all applications: from quality control of incoming goods and final products to at-line quality assurance of intermediate or final products.

- Ergonomic design and light weight for analysis on the go
- Measure comfortably using the joystick to control the automatic sampling pump, and avoid the added strain of manual sampling
- Set up a range of acceptable limits for your products
- Comprehensive result storage
- Guided workflows
- Measurement scales: API, baumé, brix, density, ethanol (alcohol), H_2SO_4 , proof (US and IP), plato, specific gravity, and user defined concentration
- The DensitoPro comes with an RFID and barcode reader, allowing for a simple point-and-shoot sample ID entry

Specifications	
Measurement Range, g/cm^3	0 to 3
Accuracy, $\pm \text{g/cm}^3$	0.001
Repeatability, g/cm^3	0.0005
Resolution, g/cm^3	0.0001
Sample Temp Range, $^\circ\text{C}$	0 to 50
Min Sample Volume, mL	2

Code	Alt Ref	Description	Price
DEN1014	30330857	Densito	£2470.00
DEN1016	30330858	DensitoPro	£2765.00

Accessories

Code	Alt Ref	Description	Pack	Price
DEN1018	30451628	EasyDirect density and refractometry PC software	Each	£615.00
DEN1020	30330860	Protective cover	Each	£120.50
DEN1022	51338010	Water density/refractive index combined standard	9mL	£295.00
DEN1024	51338012	Dodecane density/refractive index combined standard	9mL	£181.00
DEN1026	51338013	2-4-dichlorotoluene density/refractive index combined standard	9mL	£179.00
DEN1028	30449268	RFID SmartSample tags	10	£39.30
DEN1030	30449269	RFID SmartSample tag film for tanks	10	£47.20
DEN1032	30330861	Carrying case	Each	£142.50
DEN1034	30330864	Bluetooth printer Godex MX20 with dongle	Each	£605.00
DEN1036	30330863	Bluetooth dongle	Each	£142.50
DEN1038	30330865	Bluetooth printer paper	5	£11.70
DEN1040	51191926	USB cable (1.8m A-B) replacement cable to connect Density2Go to printer	Each	£18.50
DEN1042	30330847	Filling tube, 190mm	Each	£35.70
DEN1044	30330848	Filling tube, 600mm	Each	£54.00