

## Corning® CoolRacks™ Continued



Code	Description	Capacity	Colour	Price
<b>432064</b>	CoolRack 250mL	1 x 250mL Easy Grip storage bottle	Grey	<b>£335.00</b>
<b>432065</b>	CoolRack VS13	9 x 13 x 75mm blood tubes	Grey	<b>£309.00</b>
<b>432066</b>	CoolRack V16	9 x 16 x 100mm blood tubes	Grey	<b>£309.00</b>
<b>432067</b>	CoolRack V13	9 x 13 x 100mm blood tubes	Grey	<b>£309.00</b>
<b>432069</b>	CoolRack LV lightweight insulated module	12 x 13/16mm diameter tubes	Grey	<b>£292.00</b>
<b>432068</b>	CoolRack L lightweight insulated module	12 x 15mL centrifuge tubes	Grey	<b>£292.00</b>

## Corning® CoolBox™ Systems



Corning® CoolBox™ Systems are bench top cooling systems that provide sample cooling or freezing without ice, electricity, or batteries. Versatile and accommodating to a variety of sample formats and temperatures, the CoolBox™ modular design enables the use of both Corning® CoolRack® and Corning® CoolSink® thermo-conductive sample modules to hold microcentrifuge tubes, cryogenic vials, PCR tubes, or assay plates at a stable temperature while you work.

Corning® XT Starter is a low profile, open-platform cooler that accommodates most Corning® CoolRack® and Corning® CoolSink® modules, making it a versatile and flexible tool for a variety of applications. XT Starter keeps samples cold (0.5 to 4.0°C) for up to four hours. The 1 to 8°C temperature indicator provides visual assurance of proper temperature.

Corning® CoolBox™ XT systems are ideal as a benchtop cooling workstation. This system includes a base, collar, lid, and one or two Corning® XT Cooling Cores which can be used for sample cooling and accommodates one tube or plate sample module.



Code	Description	Colour	Supplied with	Price
<b>432014</b>	XT starter, complete	Purple	XT holder, XT cooling core	<b>£221.00</b>
<b>432015</b>	XT starter, holder only	Purple	XT holder	<b>£68.00</b>
<b>432021</b>	CoolBox XT	Purple	Holder, cooling core	<b>£537.00</b>
<b>432022</b>	CoolBox XT	Green	Holder, cooling core	<b>£537.00</b>
<b>432023</b>	CoolBox XT	Orange	Holder, cooling core	<b>£537.00</b>
<b>432024</b>	CoolBox XT	Pink	Holder, cooling core	<b>£537.00</b>
<b>432025</b>	CoolBox XT	Purple	Holder, 2 x cooling cores	<b>£825.00</b>
<b>432026</b>	CoolBox XT	Green	Holder, 2 x cooling cores	<b>£825.00</b>
<b>432027</b>	CoolBox XT	Orange	Holder, 2 x cooling cores	<b>£825.00</b>
<b>432028</b>	CoolBox XT	Pink	Holder, 2 x cooling cores	<b>£825.00</b>

## Accessories

Code	Description	Colour	Price
<b>432081</b>	XT cooling core for CoolBox XT Workstations for 0.5 to 4°C	—	<b>£178.00</b>
<b>432082</b>	XT freezing core for CoolBox XT Workstations for below 0°C	—	<b>£244.00</b>
<b>432083</b>	Extension collar for CoolBox XT	Purple	<b>£94.00</b>
<b>432084</b>	Extension collar for CoolBox XT	Green	<b>£94.00</b>
<b>432085</b>	Extension collar for CoolBox XT	Orange	<b>£94.00</b>
<b>432086</b>	Extension collar for CoolBox XT	Pink	<b>£94.00</b>
<b>432088</b>	Extension collar for CoolBox 2XT	Green	<b>£117.00</b>
<b>432089</b>	Extension collar for CoolBox 2XT	Orange	<b>£117.00</b>
<b>432090</b>	Extension collar for CoolBox 2XT	Pink	<b>£117.00</b>
<b>432091</b>	1 to 8°C adhesive strip LCD temperature display	—	<b>£50.00</b>
<b>432092</b>	Elastic sleeves for CoolRack modules	—	<b>£28.00</b>