

LMS Series 2 Cooled Incubators Continued

LMS

Single Set Point Control

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5324	210	135	3	600 x 660 x 945	£3169.00
INC5326	220	200	4	600 x 640 x 1370	£3311.00
INC5328	230	290	4	600 x 650 x 1800	£3580.00
INC5330	240	420	5	750 x 780 x 1770	£4402.00

Automatic Dual Temperature Control

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5332	210A	135	3	600 x 660 x 945	£3579.00
INC5334	220A	200	4	600 x 640 x 1370	£3721.00
INC5336	230A	290	4	600 x 650 x 1800	£3990.00
INC5338	240A	420	5	750 x 780 x 1770	£4812.00

Fully Programmable Controller (Multiple Temperatures/Time Frames)

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5340	210NP	135	3	600 x 660 x 945	£3519.00
INC5342	220NP	200	4	600 x 640 x 1370	£3661.00
INC5344	230NP	290	4	600 x 650 x 1800	£3930.00
INC5346	240NP	420	5	750 x 780 x 1770	£4752.00



LMS Series 3 Cooled Incubators

LMS

Offering reliable and accurate temperature-control, operating between -10 to +50°C. The range comprises of 4 cabinet sizes to suit most workloads, ranging from 90 to 450L capacity. Fitted with a stainless steel interior.

- Digital set/display
- Fan-assisted circulation
- Hermetically-sealed refrigeration system
- Self-evaporating condensation drip tray
- Preset overheat protection
- For 220/240V, single phase supply
- Stainless steel interior
- A full range of upgrades including fruit fly modification are also available to suit your requirements, please enquire for details

Single Set Point Control

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5008	100	90	2	610 x 630 x 800	£3144.00
INC5010	R3/200W	227	3	680 x 670 x 1550	£4672.00
INC5012	R3/300W	290	3	680 x 670 x 1810	£5063.00
INC5014	400	450	3	690 x 920 x 1810	£4801.00