

## LMS Series 3 Cooled Incubators Continued

**LMS**

### Automatic Dual Temperature Control

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5090	100A	90	2	610 x 630 x 800	£3554.00
INC5092	200A	227	3	680 x 670 x 1550	£4672.00
INC5094	300A	290	3	680 x 670 x 1810	£5063.00
INC5096	400A	450	3	690 x 920 x 1810	£5211.00

### Fully Programmable Controller (Multiple Temperatures/Time Frames)

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5009	100NP	90	2	610 x 630 x 800	£3494.00
INC5011	200NP	227	3	680 x 670 x 1550	£4612.00
INC5013	300NP	290	3	680 x 670 x 1810	£5003.00
INC5015	400NP	450	3	690 x 920 x 1810	£5151.00

## LMS Series 4 Cooled Incubators

**LMS**

Offering reliable and accurate temperature control, operating between -10 to +50°C. The design of these cabinets facilitates the addition of optional extras to provide even greater flexibility for individual applications. The range comprise of 2 cabinet sizes, 600 and 1200L capacity. These cabinets offer the ability to cope with large workloads and the flexibility with the addition of optional extras to suit most applications.

- Digital set/display
- Fan-assisted circulation
- Hermetically-sealed refrigeration system
- Self-evaporating condensation drip tray
- Preset overheat protection
- For 220/240V, single phase supply
- Stainless steel interior and exterior
- A full range of upgrades including fruit fly modification are also available to suit your requirements, please enquire for details



### Single Set Point Control

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5020	R4/600	600	5	800 x 900 x 2080	£5995.00
INC5022	R4/1200	1200	10	1500 x 900 x 2080	£7923.00

### Automatic Dual Temperature Control

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5100	600A	600	5	800 x 900 x 2080	£6405.00
INC5102	1200A	1200	10	1500 x 900 x 2080	£8333.00

### Fully Programmable Controller (Multiple Temperatures/Time Frames)

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5021	600NP	600	5	800 x 900 x 2080	£6345.00
INC5023	1200NP	1200	10	1500 x 900 x 2080	£8273.00