

A&D GX-A/GF-A Series Precision Balances Continued



Code	Alt Ref	Capacity, g	Readability, g	Pan Dims, mm	Calibration	Price
BAL1388	GX403A	420	0.001	128 x 128	Internal	£1597.00
BAL1394	GF603A	620	0.001	128 x 128	External	£1488.00
BAL1392	GX603A	620	0.001	128 x 128	Internal	£1854.00
BAL1398	GF1003A	1100	0.001	128 x 128	External	£1569.00
BAL1396	GX1003A	1100	0.001	128 x 128	Internal	£1957.00
BAL1356	GF1202A	1200	0.01	165 x 165	External	£1028.00
BAL1402	GF1603A	1620	0.001	128 x 128	External	£1677.00
BAL1400	GX1603A	1620	0.001	128 x 128	Internal	£2112.00
BAL1360	GF2002A	2200	0.01	165 x 165	External	£1136.00
BAL1358	GX2002A	2200	0.01	165 x 165	Internal	£1442.00
BAL1364	GF3002A	3200	0.01	165 x 165	External	£1217.00
BAL1362	GX3002A	3200	0.01	165 x 165	Internal	£1545.00
BAL1368	GF4002A	4200	0.01	165 x 165	External	£1271.00
BAL1366	GX4002A	4200	0.01	165 x 165	Internal	£1597.00
BAL1350	GF6001A	6200	0.1	165 x 165	External	£1185.00
BAL1348	GX6001A	6200	0.1	165 x 165	Internal	£1288.00
BAL1372	GF6002A	6200	0.01	165 x 165	External	£1509.00
BAL1370	GX6002A	6200	0.01	165 x 165	Internal	£1854.00
BAL1354	GF10001A	10200	0.1	165 x 165	External	£1293.00
BAL1352	GX10001A	10200	0.1	165 x 165	Internal	£1437.00
BAL1374	GX10002A	10200	0.01	165 x 165	Internal	£1906.00
BAL1376	GF10002A	10200	0.01	165 x 165	External	£1725.00

OHAUS Scout® Portable Balances



Ready for your lab, wherever that may be. With the updated Scout® portable balances, you get the same exceptional accuracy and reliability of earlier Scout® balances along with several new improvements. Easy to use with a 4-button keypad and intuitive, visual menu prompts.

- Applications include weighing (7 units), counting, percent, display hold, totalisation, check weigh, specific gravity/density, dynamic
- LCD screen with white LED backlight
- Mains or battery
- Integrated overload protection system
- ABS plastic housing with stainless steel pan



Code	Alt Ref	Model	Capacity, g	Readability, g	Pan Dims, mm	Price
BAL0030	30253070	SKX123	120	0.001	Ø 93	£408.00
BAL0031	30253071	SKX222	220	0.01	Ø 120	£272.00
BAL0032	30253072	SKX422	420	0.01	Ø 120	£330.00
BAL0033	30253073	SKX622	620	0.01	Ø 120	£363.00
BAL0034	30253074	SKX1202	1200	0.01	170 x 140	£414.00
BAL0035	30253075	SKX2202	2200	0.01	170 x 140	£481.00
BAL0036	30253076	SKX421	420	0.1	Ø 120	£193.50
BAL0037	30253077	SKX621	620	0.1	170 x 140	£222.00
BAL0038	30253078	SKX2201	2200	0.1	170 x 140	£290.00
BAL0039	30253079	SKX6201	6200	0.1	170 x 140	£369.00
BAL0040	30253080	SKX8200	8200	1	170 x 140	£311.00

Accessories

Code	Alt Ref	Description	Pack	Price
BAL1017	30268982	RS232 interface kit	Each	£55.70
BAL1019	30268987	Stacking kit	6	£50.10
BAL1018	30268985	Bluetooth interface kit	Each	£104.50