


Heathrow Scientific Smartpicks® Inoculation Picks



- The 'smart' way to isolate individual colonies
- Naturally biodegradable birch wood
- Autoclavable
- Available in semi or needle point styles for delicate operations
- The mini fits neatly into a 1.5mL microfuge tube allowing the lid to be closed for full sample transfer/recovery

Code	Alt Ref	Description	Length, mm	Pack	Price
IN00002	HS120010	Mini needle point	32.5	2000	£147.00
IN00003	HS120011	Standard semi point	63.5	1500	£162.00
IN00004	HS120012	Standard needle point	63.5	1500	£162.00
IN00006	HS120014	Large needle point	152.4	500	£90.00



Artery Forceps, Halstead's

Stainless steel Halstead's mosquito artery clamps, all 125mm long.

Code	Description	Price
INS4220	Straight clamp with box joint	£8.80
INS4222	Straight clamp with screw joint	£10.88
INS4224	Curved clamp with box joint	£12.64
INS4226	Curved clamp with screw joint	£12.38



Artery Forceps, Spencer Wells

Straight, stainless steel, Spencer Wells artery clamps.

Code	Description	Length, mm	Price
INS4240	Screw joint	130	£9.67
INS4244	Screw joint	150	£9.88
INS4248	Screw joint	175	£18.48
INS4252	Screw joint	200	£22.60
INS4242	Box joint	130	£9.40
INS4243	Box joint	150	£10.70
INS4250	Box joint	200	£18.40



Thermo Scientific™ Nalgene™ Scissor-Type Forceps



Chemically compatible with most lab reagents. Double-action hinges provide extra leverage and long life. Feature a scissor-type handle with a ratchet to lock tips shut and a handy clamp for chromatography columns. Polypropylene. Autoclavable.

Code	Tips	Pack	Price
6320-0010	Straight, serrated	12	£32.30

