



## Stainless Steel Scalpels

| Code    | Length, mm | Price |
|---------|------------|-------|
| INS4580 | 38         | £8.20 |
| INS4582 | 43         | £8.20 |
| INS4584 | 52         | £8.20 |



## Swann-Morton Disposable Scalpels

Swann-Morton LIMITED

Sterile, disposable scalpels, consisting of stainless steel blades fitted to polystyrene handles, offering a convenient alternative to the traditional blade fitted to a metal handles.

- Graduated measure on the reverse of the handle (for indication only)
- Marked with the CE symbol
- Individually packaged into a transparent thermos-formed tray
- When peeling the packet the seal changes colour to assure the user of sterility up to the point of usage

Supplied in packs of 10



| Code    | Alt Ref | Blade Pattern No. | Price |
|---------|---------|-------------------|-------|
| INS4761 | 0501    | 10                | £7.04 |
| INS4762 | 0502    | 10A               | £7.04 |
| INS4759 | 0503    | 11                | £7.04 |
| INS4757 | 0505    | 15                | £7.04 |
| INS4763 | 0506    | 20                | £7.56 |
| INS4765 | 0507    | 21                | £7.56 |
| INS4767 | 0508    | 22                | £7.56 |
| INS4755 | 0509    | 22A               | £7.56 |
| INS4753 | 0510    | 23                | £7.56 |
| INS4769 | 0511    | 24                | £7.56 |

## Swann-Morton No. 3 and No. 7 Ranges

Swann-Morton LIMITED

### Handles

For use with No. 6, 9, 10, 10A, E/11, 11, 11P, 12, 12D, 13, 14, 15, 15A, 15C, 15T, 16, 40 and SG3 blades.



| Code    | Alt Ref | Description | Length, mm | Graduated | Price |
|---------|---------|-------------|------------|-----------|-------|
| INS4603 | 0933    | No. 3       | 120        | Yes       | £6.09 |
| INS4605 | 0907    | No. 7       | 160        | No        | £7.14 |