

## Swann-Morton Qlicksmart Scalpel Blade Removers

### BladeFLASK Scalpel Blade Remover

BladeFLASK provides a simply safer way to remove and contain used scalpel blades and prevent the dangerous injuries commonly caused by removing blades with fingers, forceps or re-sheathing. A puncture proof sharps container which, when it is full, can be disposed of safely without any further exposure to sharps contact. A built-in counter with shut-off mechanism activates automatically after 100 blades have been removed.

- No danger of overfilling
- Ergonomic simplicity
- Requires safe, single-handed technique
- Compatible scalpel blades: 10, 10A, 11, 12 etc. up to and including 25A, 27 and 36
- Compatible scalpel handles: 3, 3L, 4, 4L, 5, 7 and 9
- Can be wall or bench mounted



Code	Alt Ref	Capacity	Price
INS4698	QFYUKSM	100 blades	£30.70

### BladeFLASK EVO Scalpel Blade Remover

Able to accommodate even more types of BP scalpel handles and blades than any other remover, giving you greater freedom in using any scalpel you want. It is so simple to use with just one hand. Embrace safety with the new standard in scalpel blade removals.

- Single-handed blade removal
- Shuts off after maximum capacity is reached
- Puncture-proof
- Tested to meet UN Transport Requirements for Dangerous Goods
- Compatible with all scalpel blades, including PM60A, PM60B ,and PM70
- Compatible with most standard BP scalpel handles including barron, hexagonal, bulbous, ergonomic, PM60, and PM8 handles
- Can be wall or bench mounted



Code	Alt Ref	Capacity	Price
INS4722	QBVYUKSM	100 blades	£53.00

### BladeFLASK Mounting Bracket

Universal mounting bracket allows for single-handed blade removal operation.

- Can be mounted on a wall, or on a trolley or bench
- Affixed using screws or extra strong double-sided adhesive

Code	Alt Ref	Price
INS4697	QBRGEN	£10.29

