



SLS Select Sweetie Jars



- With blue lid

Code	Capacity, L	Pack	Price
SLS7800	3	35	£86.00
SLS7802	4.5	27	£87.00



Thermo Scientific™ Nalgene™ Wide Mouth Jars



- Wide mouth with straight sides to provide full-access and complete sample retrieval
- Clear PMP with white polypropylene screw closure
- Lighter in weight than glass
- Lose only 2% of formaldehyde per year
- Recessed tops for stacking
- Good clarity with excellent chemical resistance
- Note: not liquid leakproof, not recommended for liquid transport

Supplied in packs of 4

Code	Alt Ref	Capacity, mL	Closure Size, mm	Dims, Ø x h, mm	Price
JAR2104	2117-0060	60	53	48 x 45	£17.90
JAR2106	2117-0125	125	70	64 x 74	£30.50
JAR2108	2117-0250	250	70	75 x 119	£35.60
JAR2112	2117-1000	1000	120	112 x 151	£69.00



Snapsecure Jars

- Polypropylene
- Wide mouth with straight sides
- Press-on lid

Code	Capacity, mL	Dims, Ø x h, mm	Pack	Price
JAR3002	100	49 x 58	360	£211.00
JAR3004	160	49 x 92	220	£199.00
JAR3008	265	62 x 90	144	£188.00
JAR3013	330	62 x 113	96	£157.00
JAR3014	450	75 x 106	105	£181.00
JAR3016	530	75 x 130	70	£134.00
JAR3018	650	75 x 156	70	£124.00