


Masking Tape

- General purpose paper masking tape
- Available in 2 sizes

Code	Dims, mm x m	Price
TAP1100	25 x 50	£3.00
TAP1105	50 x 50	£5.24



Tough-Spots® Cryogenic Labels



- PVC
- Accept writing from most lab markers
- Circular
- For storage in temperatures from -196 to 80°C
- Please contact us for the full range of colours and formats
- Also available in printable sheets



Supplied in rolls of 1000 labels

Code	Size, mm	For use with	Colour	Price
MIC5054	9,5	1.5mL tubes/vials	White	£39.50
MIC5052	13	2.0mL tubes/vials	White	£39.50
MIC5014	13	2.0mL tubes/vials	Blue	£40.60
MIC5016	13	2.0mL tubes/vials	Green	£49.10
MIC5042	13	2.0mL tubes/vials	Orange	£37.80
MIC5018	13	2.0mL tubes/vials	Red	£39.50
MIC5020	13	2.0mL tubes/vials	Yellow	£39.50

Tough-Tags® Labels



- Chemically inert labels that strongly adhere to microtubes, test tubes, storage boxes and other microscope slides
- Temperature resistant from -40 to 121°C
- Withstand autoclaving, boiling water baths, freezer temperatures, organic solvents and caustic agents
- For best adhesion it is recommended that siliconised surfaces be avoided
- Perfect for accurately identifying important samples stored in tubes, microplates or boxes that will be autoclaved or stored in moderate freezer temperatures
- Please contact us for the full range of colours and formats
- Also available in printable sheets



Supplied in rolls of 1000 labels

Code	Size, mm	For use with	Colour	Price
MIC5078	24 x 13	0.5mL tubes	Clear	£52.00
MIC5082	24 x 13	0.5mL tubes	Green	£52.00
MIC5084	24 x 13	0.5mL tubes	White	£62.00
MIC5086	24 x 13	0.5mL tubes	Yellow	£52.00
MIC5088	33 x 13	1.5 - 2mL tubes	Blue	£38.90
MIC5092	33 x 13	1.5 - 2mL tubes	Green	£48.00
MIC5094	33 x 13	1.5 - 2mL tubes	Yellow	£45.00
MIC5096	33 x 13	1.5 - 2mL tubes	Red	£41.00
MIC5098	33 x 13	1.5 - 2mL tubes	White	£43.00