

OHAUS Compass™ CX Balances



Designed with energy efficiency in mind, the Compass™ CX series offers a battery life of up to 1000 continuous hours, allowing for uninterrupted operation. Featuring a slim, stackable and lightweight design, the CX fits ideally in lab, industrial, education, and even household settings. The streamlined two button operation means that the balance can be operated easily, and large backlit LCD display enables easy viewing of weighing results.

- LCD display with backlight
- Powered by 3 x AA batteries (included) or AC adaptor (not included)
- ABS housing, stainless steel pan, transportation lock, slip resistant feet
- Design features include stability indicator, mechanical and software overload/underload protection, low battery indicator, auto shut off
- Slim, space saving design allows for stackable storage for up to four scales



Code	Alt Ref	Model	Capacity, g	Readability, g	Price
BAL1454	30428208	CX221	220	0.1	£67.00
BAL1456	30428209	CX621	620	0.1	£72.00
BAL1458	30428210	CX1201	1200	0.1	£78.00
BAL1460	30428211	CX2200	2200	1	£76.00
BAL1462	30428212	CX5200	5200	1	£72.00

SLS Lab Pro Moisture Analysers



- Fast and uniform heating with halogen lamp and innovative Secondary Radiation Assist (SRA) technology
- Accuracy verification with Sodium Tartrate Dihydrate
- High repeatability with the Super Hybrid Sensor (SHS) technology
- Calibration function for heating temperature (SLS5552 and SLS5554 only)
- Five different measurement programs
- Selectable heating patterns
- Translucent window enabling progress check
- Memory function for measurement conditions
- Large bright vacuum fluorescent display
- GLP/GMP/GCP/ISO compliant output
- Quick reference card provided



Code	Capacity, g	Readability, g	Price
SLS5552	71	0.0001	£4496.00
SLS5554	51	0.001	£2238.00
SLS5556	51	0.002	£2168.00

Accessories

Code	Description	Price
BAL2800	Multifunctional compact printer	£649.00