

## Cryo-Babies® and Cryo-Tags® Cryogenic Labels



- Withstand both long term cryogenic storage conditions and conventional freezer storage
- Chemically inert with a solvent resistant adhesive that will adhere to cryogenic vials, microtubes, test tubes, glass vials, beakers and metal racks
- White
- Please contact us for the full range of colours and formats



### On a Roll

- For temperature ranges from -196 to 80°C
- Withstand flash/snap freezing in liquid nitrogen

Code	Description	Size, mm	For use with	Pack	Price
<b>CRY8520</b>	Cryo-Babies®	24 x 13	0.5mL tubes	Roll of 1000 labels	<b>£52.90</b>
<b>CRY8500</b>	Cryo-Babies®	33 x 13	1.5mL tubes	Roll of 1000 labels	<b>£44.70</b>
<b>CRY8502</b>	Cryo-Tags®	38 x 19	General purpose	Roll of 1000 labels	<b>£44.80</b>



### Formatted for Laser Printers

- For temperature ranges from -196 to 150°C
- Withstand flash/snap freezing in liquid nitrogen, boiling water baths at 100°C and dry heat up to 150°C

Code	Description	Size, mm	For use with	Pack	Price
<b>CRY8512</b>	Cryo-Babies®	24 x 13	0.5mL tubes	2380 (20 sheets of 119 labels)	<b>£125.00</b>
<b>CRY8510</b>	Cryo-Babies®	33 x 13	1.5 - 2.0mL tubes	1700 (20 sheets of 85 labels)	<b>£118.00</b>
<b>CRY8514</b>	Cryo-Tags®	38 x 19	General purpose	1200 (20 sheets of 60 labels)	<b>£126.00</b>
<b>CRY8536</b>	Cryo-Tags®	43 x 19	Cryogenic vials	1040 (20 sheets of 52 labels)	<b>£99.00</b>
<b>LCRY-1000</b>	Cryo-Tags®	A4	Cut to size	20 (uncut A4 sheets)	<b>£94.00</b>



## Heathrow Scientific Label Station



- Organise your labels for easy accessibility
- Portable label station holds up to six standard size boxes of Tough-Tags®, Cryo-Babies® and Tough-Spots®
- Just drop the box in the label station and feed the paper through the slot
- Writing surface
- Cutting edge for easy dispensing
- Non-skid rubber feet
- Counter-sunk holes provided for wall mounting
- Labels not included

Code	Alt Ref	Dims, w x d x h, mm	Colour	Price
<b>LAB0010</b>	HS23460	280 x 120 x 89	Translucent blue	<b>£94.00</b>