

Neurothesiometer Continued

Accessories


Code	Description	Price
NEU1520	Spare battery pack	£122.00
NEU1522	Battery charger, UK	£131.00

McMaster Egg Slide Counting Chamber

Used for determining faecal egg counts, which is a useful indicator of an animal's health.

- Each of the two chambers has a 10 x 10mm grid permanently etched on the underside of the top cover
- Volume under each grid: 0.15mL
- Cell depth: 1.5mm

Code	Dims, w x l, mm	Price
HAE2030	76 x 32	£152.00



Cole-Parmer® Stuart™ SMP10 Digital Melting Point Apparatus

Designed for maximum safety and ease of use. Select, measure, and display the temperature digitally, ensuring accuracy and avoiding the need for a glass thermometer. Test two samples simultaneously, viewing the melt via a magnifying lens with built in illumination. Set a plateau temperature for initial warm up – unattended warm up time frees up the user until the plateau temperature is reached. Heating rate continues at 2°C after the warm up, slowly reaching the melt temperature for greater accuracy.

- Easily view the melt via magnifying lens with built in illumination
- Be more productive – test two samples simultaneously
- Extend the back feet for the optimal viewing angle
- 3 digit LED display
- Supplied with 100 x capillary tubes (open at one end)



Code	Alt Ref	No. of Samples	Temp Range, °C	Heat Up Time, °C	Dims, w x d x h, mm	Price
MEL2104	03011-41	2	Ambient to 300	20/minute to plateau, 2/minute to melt	160 x 220 x 170	£732.00

A world of products in one place...

- Over 500,000 products
- Easy ordering capability
- Detailed product information
- Promotions and offers

