



Thermo Scientific™ Oxoid™ AnaeroBox™

Thermo
SCIENTIFIC

The AnaeroBox™ is lightweight yet robust with an uncomplicated design that makes it easy to maintain and load with a compartment for both the gas kit and also indicator, minimising the possibility of malfunction. The shape also makes it easy to store and will save space in the incubator over more conventional round jars with bulky lids that include valves and gauges. The practicality of the AnaeroBox™ complements the innovative technology which makes AnaeroGen™ and CampyGen™ sachets so easy to use. The polycarbonate lid is secured to the base by four clips that make it easy to open and close.

Code	Volume, L	Capacity	Price
AB0025A	2.5	12 x 90mm or 21 x 55mm plates	£287.00

Accessories

Code	Alt Ref	Description	Pack	Price
ANA1044	AN0025A	AnaeroGen™ anaerobic gas generating sachets	10 x 2.5L sachets	£44.00
ANA1037	CN0025A	CampyGen™ gaseous atmosphere generating sachets	10 x 2.5L sachets	£32.00
CD025A	CD0025A	CO ₂ Gen™ carbon dioxide atmosphere generating sachets	10 x 2.5L sachets	£33.00
ANA1016	BR0055B	Resazurin anaerobic indicator strips	100 strips	£75.00



Thermo Scientific™ Oxoid™ AnaeroJar™

Thermo
SCIENTIFIC

Atmospheric generation system. The 2.5L jar will hold up to 12 disposable plastic Petri dishes. The jar's polycarbonate base is secured to the lid by four spring loaded clips. The clips themselves allow self-venting in the unlikely event of a pressure build up. A carrying handle is also provided, enabling the jar to be transported safely between the bench and the incubator. An easy to operate vacuum relief screw permits the inlet of air to the jar, thus overcoming any difficulty which could be caused by the presence of negative pressure. The AnaeroJar™ is lightweight and robust with excellent space saving features; there are no external clamps or protruding valves and the carrying handle tucks away into the lid of the jar itself.

Code	Alt Ref	Volume, L	Capacity	Price
ANA1040	AG0025A	2.5	12 x Petri dishes	£348.00

Accessories

A wide range of Oxoid™ and Difco™ media are available, full catalogues available on request.

Code	Alt Ref	Description	Pack	Price
ANA1044	AN0025A	AnaeroGen™ anaerobic gas generating sachets	10 x 2.5L sachets	£44.00
ANA1037	CN0025A	CampyGen™ gaseous atmosphere generating sachets	10 x 2.5L sachets	£32.00
ANA1050	AG0026A	AnaeroJar™ base	Each	£162.00
ANA1056	AG0029A	AnaeroJar™ plate carrier	Each	£186.00
ANA1058	AG0030A	AnaeroJar™ 'O' ring	5	£42.00
AN0020C	AN0020C	Replacement sachets for AnaeroGen™ Compact gas generation system	20	£74.00
CD025A	CD0025A	CO ₂ Gen™ carbon dioxide atmosphere generating sachets	10 x 2.5L sachets	£33.00
AN010C	AN0010C	AnaeroGen™ compact, 10 sachets and 10 pouches	10	£61.00