


E&O Laboratories Animal Blood Products



Freshly harvested aseptically collected donor horse blood and sheep blood. These blood products are oxygenated and supplied in both bottles and bags for convenient aseptical addition to culture media. Plasma products can be fractionated after the addition of a range of anticoagulants: citrate, Alsever's, sodium oxalate, heparin and EDTA. All donor animal blood products are harvested from stress free animals with guaranteed PCVs and free from any microbiological inhibitors. The most common use is for the manufacture of blood agar plates.



Code	Alt Ref	Description	Packaging	Pack	Price
MED1300	DHB025-P030-25	Defibrinated horse blood	Universal polystyrene bottle	25mL	£5.80
MED1302	DHB100-Y125-100	Defibrinated horse blood	Alpha syrup bottle with 28mm polypropylene cap	100mL	£12.40
MED1306	DHB500-N500-500	Defibrinated horse blood	PETG media bottle with 28mm PETG tamper evident cap	500mL	£57.00
MED1304	DHB1L-1LB-1000	Defibrinated horse blood	Media bag	1L	£95.00
MED1308	DSC025-P030-25	Defibrinated sheep blood	Universal polystyrene bottle	25mL	£8.20
MED1310	DSC100-Y125-100	Defibrinated sheep blood	Alpha syrup bottle with 28mm polypropylene cap	100mL	£21.80
MED1314	DSC500-N500-500	Defibrinated sheep blood	PETG media bottle with 28mm PETG tamper evident cap	500mL	£90.00
MED1312	DSC1L-1LB-1000	Defibrinated sheep blood	Media bag	1L	£147.00
MED1316	HBL100-Y125-100	Horse blood lysed	Alpha syrup bottle with 28mm polypropylene cap	100mL	£17.10
MED1318	HIS100-N100-100	Heat inactivated horse serum	PETG media bottle with 28mm PETG tamper evident cap	100mL	£17.00
MED1320	SCA025-P030-25	Sheep cells in alsevers	Universal polystyrene bottle white polystyrene cap	25mL	£7.47
MED1322	SHS100-N100-100	Sterile horse serum	PETG media bottle with 28mm PETG tamper evident cap	100mL	£13.40



BD Culture Media and Ingredients



Agar

Code	Description	Brand	Pack	Price
212304	Grade A	BBL™	500g	£162.00
214530	Bacteriological	Difco™	500g	£141.00
214510	Bacteriological	Difco™	2kg	£423.00
214520	Bacteriological	Difco™	10kg	£1893.00
214220	Noble	Difco™	100g	£205.00
214230	Noble	Difco™	500g	£417.00
214050	Purified	(Bacto™) Difco™	100g	£95.00
214010	Purified	(Bacto™) Difco™	454g	£192.00
214030	Purified	(Bacto™) Difco™	2kg	£751.00
214040	Purified	(Bacto™) Difco™	10kg	£3135.00
281230	Technical	Difco™	500g	£229.00
281210	Technical	Difco™	2kg	£501.00
214340	Gelatin	Difco™	500g	£89.00