

Gilson TRACKMAN® Duo-Plate Tracker



A highly effective simple to use tool that helps avoid costly and time consuming pipetting mistakes and prevent cross contamination. Designed to accommodate most 96 and 384 well plates, TRACKMAN® keeps track of where you are pipetting through bright red LED markers and can be programmed to accommodate up to 16 different pipetting patterns. Through the use of a microtube accessory, TRACKMAN® can also assist in the transfer of liquids between microtubes and multiwell plates.

- Accommodates two multiwell plates side by side
- Clear LED marker to indicate current well/s
- Fully programmable and can accommodate pipetting patterns in rows, columns or any other desired pattern
- Intuitive software, requires no training
- Battery powered to ensure portability - rechargeable battery life of up to 8 hours
- Notches ensure plates are firmly held in place, retractable stand allows unit to be used at a slight incline



Code	Price
F70301	£766.00

Gilson TRACKMAN® Connected



A tablet with exclusive accessories and applications that make pipetting on 96 and 384 well plates faster and more reliable, improving efficiency at the bench by tracking pipetting tasks. Designed to communicate with Gilson's PIPETMAN® M Connected via Bluetooth®, the tablet interacts in real time with the pipette and guides the user through their protocol with the PipettePilot™ app. The app displays where, when, and how much to pipette and even captures notes, pictures and monitors the protocol. It can even adjust the pipette volume automatically according to the method.

- Improves traceability and reproducibility
- Accelerates collaborative reporting
- Reports on your pipetting performance
- Guides you through your pipetting protocol
- Increases your efficiency and reduces costly errors



Code	Price
FB1021	£1586.00

Accessories

Code	Description	Price
FB02201	Tablet holder	£102.00
FB02202	Plate holder	£26.00

