

Tecan Sunrise Plate Readers Continued

TECAN.

Code	Alt Ref	Type	Supplied with	Temp Range, °C	Price
MIC9250	30133441	4 filter	405, 450, 492 & 620nm filters	Ambient	£6387.00
MIC9256	30133444/10037407	4 filter	405, 450, 492 & 620nm filters	RT to 42 (resolution 0.1)	£7663.00
MIC9252	30133442	6 filter	405, 450, 492, 540, 570 & 620nm filters	Ambient	£7334.00
MIC9258	30133445/10037407	6 filter	405, 450, 492, 540, 570 & 620nm filters	RT to 42 (resolution 0.1)	£8940.00
MIC9254	30133443	Gradient filter	400 - 700nm range filters	Ambient	£8625.00
MIC9260	30133446/10037407	Gradient filter	400 - 700nm range filters	RT to 42 (resolution 0.1)	£9957.00

Accessories

Code	Alt Ref	Description	Price
MIC9262	30026966	QC-Pac 2 quality control tool for IQ/OQ	£1040.00



Tecan Infinite 200 PRO Plate Readers

TECAN.

An easy to use multimode plate reader family that offers affordable, high performance detection solutions, referenced in more than 1800 peer reviewed publications.

- Absorbance - superior performance and free wavelength choice - from ELISAs to low sample volume quantifications
- Fluorescence top & bottom - choose between the sensitivity of filters and the flexibility of monochromators
- Te-Inject reagent injectors
- Automated pathlength correction
- Absorbance cuvette port, NanoQuant Plate
- Temperature control (RT +5 to 42°C)
- Plate shaking (linear and orbital)
- QC tools for IQ/OQ services
- Luminescence - optimised for glow, flash and multicolour applications performed in 384 well plates

Code	Alt Ref	Model	Absorbance (230-1000nm)	Fluorescence (230-850nm)	Luminescence	Price
MIC9264	30133334	Infinite M Nano	Monochromator	No	No	£11945.00
MIC9266	30133335	Infinite M Nano+	Monochromator	Monochromator	No	£26390.00
MIC9268	30133438	Infinite M Plex	Monochromator	Monochromator	Yes	£29782.00
MIC9270	30133437	Infinite Lumi	No	No	Yes	£12259.00
MIC9272	30133336	Infinite F Nano+	Filter	Filter	No	£15656.00
MIC9274	30133439	Infinite F Plex	Filter	Filter	Yes	£21320.00

Accessories

Code	Alt Ref	Description	Price
MIC9276	30015898	Cuvette option (Infinite M reader only)	£1565.00
MIC9278	30033939	NanoQuant plate for low volume DNA/RNA	£2667.00
MIC9280	30124710	1x channel Te-inject short tubing	£4123.00
MIC9282	30124711	1x channel Te-inject long tubing	£4123.00
MIC9284	30124708	2x channel Te-inject short tubing	£5327.00
MIC9286	30124709	2x channel Te-inject long tubing	£5328.00