



1

YEAR

## Tecan HydroSpeed Plate Washers



An advanced microplate washer, designed to provide optimised washing of cells, beads and ELISAs in both 96 and 384 well formats. It offers full control over critical wash parameters to help maximise the efficiency of cell assays, bead assays and ELISAs. With full on-board operation via the easy-to-use touchscreen interface, the HydroSpeed really does offer you power and control for all your cell, bead and ELISA washing applications.

- Cell Protection™ – extra gentle cell washing with tunable wash power
- 5.7" touch screen
- Reduced turbulence via online optimisation of dispense position
- Typical cell recovery of ≥85% using tunable aspiration settings
- Anti Clogging™ – reliable operation through active prevention of needle clogging
- High bead recovery rates (typically ≥90%) via tunable aspiration settings
- Automated rinsing and soaking of wash head for easy handling
- Dual magnets per well – efficient bead settling and high bead recovery rates
- Full on-board operation – easy to use touchscreen interface
- Tested reliability – high quality components ensure smooth operation
- Easy X-change wash heads for straightforward ultrasonic cleaning
- 21 CFR Part 11/98/79/EC IVD-D

Code	Alt Ref	Channel	No. of Inlets	Price
<b>MIC9342</b>	30054622/30055364	96 HT	1	£9779.00
<b>MIC9344</b>	30054622/30055365	96 HT	2	£10026.00
<b>MIC9346</b>	30054622/30055366	96 HT	3	£10586.00
<b>MIC9348</b>	30054622/30055367	96 HT	4	£11158.00
<b>MIC9350</b>	30055360/30055364	96 Indexing	1	£10445.00
<b>MIC9352</b>	30055360/30055365	96 Indexing	2	£10800.00
<b>MIC9354</b>	30055360/30055366	96 Indexing	3	£11360.00
<b>MIC9356</b>	30055360/30055367	96 Indexing	4	£11933.00
<b>MIC9358</b>	30054623/30055364	384	1	£12169.00
<b>MIC9360</b>	30054623/30055365	384	2	£12525.00
<b>MIC9362</b>	30054623/30055366	384	3	£13086.00
<b>MIC9364</b>	30054623/30055367	384	4	£13657.00

### Accessories

Code	Alt Ref	Description	Price
<b>MIC9366</b>	30060862	Option 4-inlet channels large volume lid incl. 4x10L wash bottles with lid, 1x 20L waste bottle with lid and 1x 4L foam trap without lid (order with instrument only)	£4616.00
<b>MIC9368</b>	30057606	Smart-2 MBS 96 plate carrier with 2 magnets per well	£2461.00
<b>MIC9370</b>	30057605	MBS 96 plate carrier with 1 magnet per well	£1486.00
<b>MIC9372</b>	30057607	MBS 384 plate carrier with 1 magnet per well	£3062.00
<b>MIC9374</b>	30055308	Hydrospeed vacuum filtration option (order with instrument only)	£2690.00



1

YEAR



## Asys Atlantis Microplate Washers



Microplate washer for standard ANSI SBS footprint 96 well plates with flexibility in washing options, simple operation and digitally controlled aspiration and dispensing pumps. The perfect choice to compliment any microplate reader in any applications environment.

- Configurable to many different well formats for maximum flexibility
- Optional 16 way manifold for washing 384 well plates
- Detachable plate carrier for easy cleaning
- Low noise, pressure free aspiration and dispensing pumps
- Auto rinse function prevents liquid path clogging
- Reduces aerosols due to cover and pumping technology