

Moticam X³ Plus Wi-Fi Camera

Motic®

The Moticam X³ Plus transfers live images up to 30 fps and captures images at 4.0MP; a compact Wi-Fi camera to share live images with up to 6 devices at the same time. A powerful high res resolution to turn your microscope into a Wi-Fi station. Simply log on with your tablet or smartphone and get ready to capture, edit and save your images with the free MotiConnect App or use your computer with Motic Images Plus 3.0 software for further image processing.

- CMOS sensor
- Wi-Fi and ethernet (RJ45)
- 12mm focusable lens
- Use Motic Images Plus 3.0 for PC, Mac and Linux or free MotiConnect app



Code	Alt Ref	Sensor Size	Imaging Area, mm	Capture Resolution, MP (pixels)	Price
MIC6505	1100600101161	1/3"	5.44 x 3.07	4 (2592 x 1520)	£448.00

Moticam BTX8/10 Plus Tablet Cameras

Motic®

The Moticam BTX Plus combines a 4MP Wi-Fi camera with an 8/10" LCD Android tablet to deliver fast live images (up to 30 fps) on any trinocular microscope. Detaching the tablet for an even more flexible handling is possible thanks to Wi-Fi image data transfer. Use the pre-installed MotiConnect App and get ready to capture, edit and share your microscope images from the tablet in a few taps. It also comes with an ethernet port, allowing you to access it from your computer and the Motic Images Plus 3.0 software (not included).

- Still image and video capture, live and still image measurement, image adjustments, white balance, automatic and manual exposure, individual objective calibration system
- Package includes: 16GB MicroSD card, USB 2.0 cable, 4-dot calibration slide and universal power supply



Code	Alt Ref	Description	Price
MIC4170	1100600101141	Moticam BTX8 with 8" tablet	£654.00
MIC4172	1100600101151	Moticam BTX10 with 10" tablet	£937.00

Sustainable laboratory solutions!

Find out how we can help you with our sustainable laboratory solutions.

To request/download a brochure or to find out more, visit:

scientificlabs.co.uk/sustainability (UK & Northern Ireland)

scientificlabs.ie/sustainability (Republic of Ireland)

scientificlabs.com/sustainability (International)

