

Schott KL 1600 Illuminator



A standard LED cold light source for fibre optics, the KL 1600 LED is an ideal option for all objects and applications. Its multiple chip white LED light engine offers an attractive alternative to conventional halogen cold light sources.

- Patented white light engine
- Easy to use fixture mechanism with clamping screw
- Two position filter slider
- Maximum light intensity of 680lm
- Continuous dimming from 0 - 100%
- Maintenance free LED engine with 50000 operating hours



Code	Alt Ref	Luminous Flux, lm	Colour Temp, K	Cooling	Dims, w x d x h, mm	Weight, kg	Price
MIC2546	150600	680	5600	Ultra low noise fan	114 x 231 x 137	2.45	£1023.00

Accessories

Code	Alt Ref	Description	Ø, mm	Length, mm	Neck Ø, mm	Price
MIC2550	400.053	UK power cord, 1.8m	—	—	—	£13.01
MIC2552	154.101	Single gooseneck light guide	4.5	600	4.5	£223.00
MIC1502	9.705152	Double gooseneck light guide	4.5	600	4.5	£332.00
MIC2566	155.204	Double branch flexible light guide	4.5	1000	4.5	£491.00

SLS Select Economy Microscope Slides



- Single frost
- 45° ground edges
- Dims, w x l: 76 x 26mm
- Thickness: 1 - 1.2mm

Code	Dims, w x l, mm	Thickness, mm	Pack	Price
SLS30500	76 x 26	1 - 1.2	50	£2.30



SLS Select Standard Microscope Slides



- Ground edge
- Available plain or twin frosted
- Dims, w x l: 76 x 26mm
- Thickness: 1mm

Supplied in packs of 50

Code	Description	Price
MIC2000	Plain	£4.00
MIC2002	Twin frost	£2.76

