


SLS Select Microscope Cover Slips



Thickness	
No. 0	0.085 - 0.130mm
No. 1	0.130 - 0.155mm
No. 1½	0.155 - 0.185mm
No. 2	0.190 - 0.250mm



Square

Supplied in packs of 200

Code	Dims, w x l, mm	Thickness	Price
MIC3100	18 x 18	No. 0	£6.15
MIC3102	20 x 20	No. 0	£7.66
MIC3104	22 x 22	No. 0	£8.29
MIC3110	18 x 18	No. 1	£5.88
MIC3112	20 x 20	No. 1	£7.66
MIC3114	22 x 22	No. 1	£8.57
MIC3116	24 x 24	No. 1	£10.19
MIC3120	18 x 18	No. 1½	£7.66
MIC3122	20 x 20	No. 1½	£9.15
MIC3124	22 x 22	No. 1½	£8.71
MIC3126	24 x 24	No. 1½	£12.54
MIC3130	18 x 18	No. 2	£8.05
MIC3132	20 x 20	No. 2	£9.15
MIC3134	22 x 22	No. 2	£8.71



Rectangular

Supplied in packs of 100

Code	Dims, w x l, mm	Thickness	Price
MIC3200	22 x 26	No. 0	£6.46
MIC3202	22 x 32	No. 0	£7.42
MIC3204	22 x 40	No. 0	£8.71
MIC3206	22 x 50	No. 0	£10.25
MIC3208	22 x 64	No. 0	£13.56
MIC3220	22 x 26	No. 1	£5.11
MIC3222	22 x 32	No. 1	£6.15
MIC3224	22 x 40	No. 1	£7.80
MIC3226	22 x 50	No. 1	£9.62
MIC3228	22 x 64	No. 1	£11.58
MIC3230	24 x 32	No. 1	£5.97
MIC3232	24 x 40	No. 1	£7.42
MIC3234	24 x 50	No. 1	£12.18
MIC3240	22 x 26	No. 1½	£5.30
MIC3242	22 x 32	No. 1½	£7.42
MIC3244	22 x 40	No. 1½	£8.35
MIC3246	22 x 50	No. 1½	£11.21
MIC3250	24 x 32	No. 1½	£7.42
MIC3252	24 x 40	No. 1½	£6.41
MIC3260	22 x 26	No. 2	£6.15
MIC3262	22 x 32	No. 2	£7.42
MIC3264	22 x 40	No. 2	£9.86
MIC3268	22 x 64	No. 2	£14.73