


Microscope Slide Staining Jars



MIC6000



MIC6014



MIC6004

Code	Description	Material	Dims, mm	Price
MIC6006	Staining dish with cavity	Soda lime glass	w x d x h: 40 x 40 x 15	£6.30
MIC6000	Coplin staining jar with lid	Soda lime glass	h x Ø: 115 x 90	£10.10
MIC6002	Square jar	PMP	h x Ø: 58 x 53	£51.00
MIC6004	Hellendahl staining jar with lid	Soda lime glass	w x d x h: 58.5 x 58.8 x 100	£26.20
MIC6014	Schiefferdecker staining jar with lid	Soda lime glass	w x d x h: 70 x 90 x 50	£17.90
MIC6008	Rectangular staining jar with lid	Soda lime glass	w x d x h: 85 x 105 x 74	£20.20
MIC6012	Rack for MIC6008	Soda lime glass	w x d x h: 67 x 87 x 50	£12.80

Simport® Scientific EasyDip™ Staining Jars



- Made of acetal, which will not break
- Will resist most staining agents including alcohol and xylene (but not phenol, iodine, or ferric chloride)
- Wide stable base offers greater stability while the inside is recessed, allowing for a smaller reagent volume of only 80mL
- Easy to clean and no metals to corrode
- Ideal for special stains, frozen sections and special processes
- Available in different colors to help better identifying contents or applications
- Dims, l x w x h: 76 x 64 x 92mm
- Can be used with 12 position slide staining rack (sold separately)



Supplied in packs of 6

Code	Alt Ref	Colour	Price
HIS0273	M900-12B	Blue	£82.00
HIS0274	M900-12G	Green	£82.00
HIS0275	M900-12P	Pink	£82.00
HIS0276	M900-12W	White	£84.00
HIS0277	M900-12Y	Yellow	£84.00

Slide Staining Rack

- Will hold up to 12 microscope slides with dimensions such as 75 x 25mm [3 x 1"] or 76 x 26mm with a thickness of 1.0 and 1.2mm
- The slides fit into individual slots for free passage and rapid drainage of staining fluids



Code	Alt Ref	Capacity	Dims, w x d x h, mm	Pack	Price
HIS0278	M905-12DGY	12 slides	60 x 64 x 97	6	£90.00