



Grant SBB Aqua Plus Boiling Baths



Robust, reliable and provide continuous 100°C operation. The range consists of four analogue models to suit a range of applications and any budget.

- Adjustable energy regulator provides steady boiling
- Constant level device maintains liquid level
- Robust and reliable design to withstand everyday wear and tear
- Simple, intuitive controls with quick and easy to set temperature
- Non drip polycarbonate lid included as standard

Code	Alt Ref	Capacity, L	Ext Dims, w x d x h, mm	Int Dims, w x l, mm	Price
BAT3200	SBBAQUA5+	5	390 x 250 x 375	145 x 290	£800.00
BAT3026	SBBAQUA12+	12	390 x 385 x 440	315 x 290	£860.00
BAT3028	SBBAQUA18+	18	390 x 600 x 445	495 x 290	£920.00
BAT3034	SBBAQUA26+	26	390 x 600 x 445	495 x 290	£990.00



Nickel Electro Clifton® Rectangular Boiling Baths



Rectangular, aluminium electrically heated, multi-place baths regulate from a gentle simmer to a vigorous boil, fitted with concentric ringed lid to vary the hole size. Features an integral constant level device to maintain constant water level during boiling and supplied with a set of rings to vary the diameter of aperture to accommodate different sizes of vessels e.g. evaporating dishes and glassware, including round bottom flasks. Fitted with a sheathed immersion element and a resettable thermal cut out for safety.

- All baths have an integral constant level device to regulate the level of water
- Supplied with concentric ringed lid to vary the diameter holes
- Smooth walled tank making it easy to clean and descale
- Features a built in energy regulator, which controls the amount of boiling in a bath from a gentle simmer to a vigorous boil
- Each bath is uniquely identified with its own serial number recorded on our product database, which includes all of its manufacturing inspection data

Code	Alt Ref	Dims, l x w x d, mm	No. of Places	Price
BAT1022	1101E	321 x 219 x 90	6 (30, 42, 58 or 80mm holes)	£539.00
BAT1024	1102E	425 x 321 x 90	6 (30, 42, 58, 80 or 101mm holes)	£574.00
BAT1026	1104E	425 x 321 x 90	12 (30, 42, 58 or 80mm holes)	£678.00

Energy Regulator

Optional accessory energy regulator can be used for bath control between vigorous boil and gentle simmer.



Code	Alt Ref	Dims, w x d x h, mm	Price
BAT1003	ER-1	128 x 125 x 80	£238.00