

SLS Lab Pro Unstirred Water Baths



Provide a stable temperature environment with digital PID temperature control for precision. Controls feature auto alarm settings, illuminated on/off switch, heater and over temperature indications.

- Sensitivity: $\pm 0.2^{\circ}\text{C}$
- Uniformity: $\pm 0.1^{\circ}\text{C}$
- Temperature range: ambient $+5$ to 99°C
- Digital PID control with LED temperature display
- Over temperature alarm with heater cut off $+4^{\circ}\text{C}$ above set point
- Touch pad controls with easy wipe surface.
- Painted body and controls features antibacterial protection, hygienic coating which actively inhibits bacterial growth
- Supplied with lid, recommended for use above 60°C



Code	Capacity, L	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Drain Outlet	Price
SLS1274	2	332 x 185 x 135	300 x 150 x 65	No	£486.00
SLS1176	4	332 x 185 x 290	300 x 150 x 150	No	£527.00
SLS1178	14	361 x 332 x 290	325 x 300 x 150	No	£619.00
SLS1180	22	538 x 322 x 290	500 x 300 x 150	Yes	£691.00
SLS1181	28	538 x 332 x 290	500 x 300d x 200	Yes	£800.00

Lids

Lift off, clear lids.

Code	For use with	Price
SLS1192	4L bath	£74.00
SLS1194	14L bath	£74.00
SLS1196	22L bath	£74.00

JULABO PURA Unstirred Water Baths



Count on proven functions in the new PURA series of water baths, but also on simple, intuitive operation and high material and component quality. The result is functional and high quality water baths that facilitate regular temperature applications in the laboratory and, thanks to their durable design, can with stand permanent loads without any problems.

- Splashproof protected mains switch
- Built in dry running protection
- Integrated timer
- Non slip rubber feet
- Ergonomic recessed grips
- Working temperature range: $+18$ to 99.9°C
- Temperature stability: $\pm 0.15^{\circ}\text{C}$
- Heating capacity up to 2kW
- Bright display
- Specially designed edge surfaces, which drain condensed fluid back into the bath
- Removable platform for full immersion of the sample containers (included)
- 1PLUS warranty offers an extension to the normal warranty of up to 24 months, limited to a maximum of 10000 working hours



Code	Alt Ref	Description	Capacity, L	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Heat Output, kW	Price
BAT1032	9550404	PURA 4	0.8 - 4.5	210 x 350 x 220	120 x 270 x 140	0.8	£886.00
BAT1034	9550410	PURA 10	1.4 - 9.5	310 x 350 x 220	220 x 270 x 140	1.2	£1014.00
BAT1036	9550414	PURA 14	2 - 14	420 x 350 x 220	330 x 270 x 140	1.8	£1050.00
BAT1038	9550422	PURA 22	3.4 - 25.5	640 x 350 x 260	550 x 270 x 180	2	£1153.00
BAT1040	9550430	PURA 30	4.8 - 36	860 x 350 x 260	770 x 270 x 180	2	£1385.00